



RETRACTOR MAKER

PRODUCT BROCHURE OF

Bose Tracheal Retractor

CERTIFIED & APPROVED

✓ EN ISO 13485:2016 + A11:2021



ISO 9001:2015



EN ISO 7153-1:2016



cGMP Approved



SCCI Registered since 2023



SIMAP Approved since 2023

All procedures and processes follow FDA guidelines.



Bose Tracheal Retractor

SKU RM-BOTR-1

\$29 USD

AVAILABLE SIZES & PRICING

SIZE / OPTION	SKU / DETAILS	PRICE (USD)
Left Orientation	RM-BOTR-1-1	\$29
Right Orientation	RM-BOTR-1-2	\$29

PRODUCT DESCRIPTION

The Bose Tracheal Retractor is a specialized airway surgical retractor engineered specifically to provide reliable, atraumatic exposure of the trachea during emergent and elective tracheostomy procedures. Designed with an anatomically curved hook blade, this instrument allows surgeons to engage and elevate tracheal rings cleanly while displacing overlying pretracheal fascia and thyroid isthmus tissue.

Built for precise directional traction, this Bose Tracheal Retractor model features an overall working length of 15 cm (\$29) and is offered in selectable Left and Right structural orientations. This dedicated directional geometry provides surgical teams with targeted lateral access to accommodate varying patient anatomy and procedural approaches, while its flat ergonomic handle ensures a steady, non-slip manual grip during crucial airway cannulation.

Manufactured from Premium Grade AISI 304/316 Stainless Steel (Non-Magnetic) with a reflection-free satin matte finish, Retractor Maker supplies the Bose Tracheal Retractor non-sterile and fully autoclavable inside protective packaging, providing ENT, trauma, and critical care units with an essential, corrosion-resistant instrument ready for dependable operative use.



THANK YOU

We look forward to discussing your surgical retractor requirements.

CONTACT RETRACTOR MAKER

EMAIL

info@retractormaker.com

WHATSAPP

+92 327 8345101

WEBSITES

retractormaker.com
retractormanufacturer.com

ADDRESS

Retractor Maker, Near FFRV+F29, Roras Road, Sialkot-51310, Pakistan

FOLLOW US



[LinkedIn](#)



[Facebook](#)



[Instagram](#)



[X \(Twitter\)](#)



[YouTube](#)



[Pinterest](#)