



RETRACTOR MAKER

PRODUCT BROCHURE OF

Bauer Sigmoid Notch Intra Oral Retractor

CERTIFIED & APPROVED

✓ EN ISO 13485:2016 + A11:2021



ISO 9001:2015



EN ISO 7153-1:2016



cGMP Approved



SCCI Registered since 2023



SIMAP Approved since 2023

All procedures and processes follow FDA guidelines.



Bauer Sigmoid Notch Intra Oral Retractor

SKU RM-BSNIOR-1

\$79 USD

AVAILABLE SIZES & PRICING

SIZE / OPTION	SKU / DETAILS	PRICE (USD)
Without Light Guide (Left or Right)	RM-BSNIOR-1-1	\$79
With Fiber-Optic Light Guide (Left or Right)	RM-BSNIOR-1-2	\$149

PRODUCT DESCRIPTION

The Bauer Sigmoid Notch Intra Oral Retractor is a specialized maxillofacial surgical retractor engineered specifically to provide secure, stable exposure of the mandibular ramus during orthognathic and reconstructive surgery. Designed with an anatomical tip contoured to seat directly into the sigmoid notch of the mandible, this instrument allows surgeons to retract the masseter muscle and periosteal envelope safely away from cutting instrumentation.

Engineered for procedural versatility, this Bauer Sigmoid Notch Intra Oral Retractor is offered in dedicated Left and Right orientation notches to match patient anatomy. It features a slightly curved blade paired with a hollow hooked handle that provides ergonomic leverage and fatigue-free control. Surgeons can choose between a standard model (\$79) and an illuminated configuration (\$149) featuring an integrated fiber-optic light guide that delivers shadowless illumination into deep intraoral corridors.

Manufactured from Premium Grade AISI 304/316 Stainless Steel (Non-Magnetic) with a glare-absorbing satin finish, Retractor Maker supplies the Bauer Sigmoid Notch Intra Oral Retractor non-sterile and fully autoclavable inside protective packaging, providing maxillofacial surgeons with high durability and surgical precision.



THANK YOU

We look forward to discussing your surgical retractor requirements.

CONTACT RETRACTOR MAKER

EMAIL

info@retractormaker.com

WHATSAPP

+92 327 8345101

WEBSITES

retractormaker.com
retractormanufacturer.com

ADDRESS

Retractor Maker, Near FFRV+F29, Roras Road, Sialkot-51310, Pakistan

FOLLOW US



[LinkedIn](#)



[Facebook](#)



[Instagram](#)



[X \(Twitter\)](#)



[YouTube](#)



[Pinterest](#)