



RETRACTOR MAKER

PRODUCT BROCHURE OF

Allport Mastoid Retractor

CERTIFIED & APPROVED

✓ EN ISO 13485:2016 + A11:2021



ISO 9001:2015



EN ISO 7153-1:2016



cGMP Approved



SCCI Registered since 2023



SIMAP Approved since 2023

All procedures and processes follow FDA guidelines.



Allport Mastoid Retractor

SKU RM-AMR-1

\$59 USD

AVAILABLE SIZES & PRICING

SIZE / OPTION	SKU / DETAILS	PRICE (USD)
Allport Mastoid Retractor	10.2 cm overall length, 4×4 sharp prongs	\$59
Protective Sleeve / Box	Individual transit packaging	\$59

PRODUCT DESCRIPTION

The Allport Mastoid Retractor is an indispensable self-retaining surgical retractor engineered specifically for otologic and skull base procedures requiring steady, hands-free exposure of the temporal bone. By holding back the retroauricular flap and periosteum with sharp interlocking prongs, this instrument frees surgical staff while providing an unencumbered microsurgical corridor.

Designed for precise tensioning in confined working fields, this Allport Mastoid Retractor configuration features an overall length of 10.2 cm (approximately 4 inches) and a 4×4 sharp prong layout. The integrated mechanical thumbscrew handle allows the operating surgeon to adjust blade expansion millimeter by millimeter, ensuring controlled counter-traction that prevents tearing of delicate periosteal tissues.

Machined from Premium Grade AISI 304/316 Stainless Steel (Non-Magnetic) with a reflection-free satin matte finish, Retractor Maker manufactures the Allport Mastoid Retractor to withstand repetitive hospital sterilization cycles, giving ENT and neurosurgical teams dependable holding power and longevity.



THANK YOU

We look forward to discussing your surgical retractor requirements.

CONTACT RETRACTOR MAKER

EMAIL

info@retractor maker.com

WHATSAPP

+92 327 8345101

WEBSITES

retractor maker.com
retractor manufacturer.com

ADDRESS

Retractor Maker, Near FFRV+F29, Roras Road, Sialkot-51310, Pakistan

FOLLOW US



[LinkedIn](#)



[Facebook](#)



[Instagram](#)



[X \(Twitter\)](#)



[YouTube](#)



[Pinterest](#)