



RETRACTOR MAKER

PRODUCT BROCHURE OF

Miskimon Cerebellar Retractor

CERTIFIED & APPROVED

✓ EN ISO 13485:2016 + A11:2021



ISO 9001:2015



EN ISO 7153-1:2016



cGMP Approved



SCCI Registered since 2023



SIMAP Approved since 2023

All procedures and processes follow FDA guidelines.



Miskimon Cerebellar Retractor

SKU RM-MCR-19-1

\$139 USD

AVAILABLE SIZES & PRICING

SIZE / OPTION	SKU / DETAILS	PRICE (USD)
Miskimon Cerebellar Retractor	Dimensions / Details: 19 cm overall length, 7×7 outward sharp prongs, ratchet lock · Quantity Included: 1	\$139
Protective Sleeve / Box	Dimensions / Details: Individual transit packaging · Quantity Included: 1	\$139

PRODUCT DESCRIPTION

The Miskimon Cerebellar Retractor is a specialized neurosurgical self-retaining instrument engineered specifically to provide reliable, hands-free scalp and muscle retraction during posterior fossa and suboccipital craniotomies. Designed with heavy-duty curved arms, this retractor anchors into thick nuchal musculature and galeal margins, establishing an expansive, static operative window over the occipital bone and cerebellar hemispheres.

Built in a dedicated mechanical footprint, this Miskimon Cerebellar Retractor configuration features an overall tool length of 19 cm (\$139) tailored for specialized wound border management. Handcrafted from Premium Grade AISI 304/316 Stainless Steel (Non-Magnetic), the working ends incorporate a 7×7 outward-curving sharp prongs array engineered to securely grip dense occipital fascia without slipping under heavy counter-traction. The base mechanical ring handle is integrated with an adjustable ratchet lock gear that holds tension reliably, eliminating assistant hand fatigue throughout extensive neurosurgical resections.

Forged with an anti-glare satin matte finish, Retractor Maker manufactures the Miskimon Cerebellar Retractor to withstand repetitive hospital autoclave cycles, providing neurosurgical operating suites and spine centers worldwide with exceptional mechanical stability, corrosion resistance, and complete clinical durability.



THANK YOU

We look forward to discussing your surgical retractor requirements.

CONTACT RETRACTOR MAKER

EMAIL

info@retractormaker.com

WHATSAPP

+92 327 8345101

WEBSITES

retractormaker.com
retractormanufacturer.com

ADDRESS

Retractor Maker, Near FFRV+F29, Roras Road, Sialkot-51310, Pakistan

FOLLOW US



[LinkedIn](#)



[Facebook](#)



[Instagram](#)



[X \(Twitter\)](#)



[YouTube](#)



[Pinterest](#)