



RETRACTOR MAKER

PRODUCT BROCHURE OF

Desai Self Retaining Retractor

CERTIFIED & APPROVED

✓ EN ISO 13485:2016 + A11:2021



ISO 9001:2015



EN ISO 7153-1:2016



cGMP Approved



SCCI Registered since 2023



SIMAP Approved since 2023

All procedures and processes follow FDA guidelines.



Desai Self Retaining Retractor

SKU RM-DSRR-1

\$129 USD

AVAILABLE SIZES & PRICING

SIZE / OPTION	SKU / DETAILS	PRICE (USD)
Desai Self Retaining Retractor	7.25" (18.41 cm) length, 90° open-loop blades, ratchet lock	\$129
Protective Sleeve / Box	Individual transit packaging	\$129

PRODUCT DESCRIPTION

The Desai Self Retaining Retractor is an innovative surgical retractor engineered specifically to provide hands-free operative field exposure while allowing direct instrumentation access through its unique blade framework. Designed with an advanced dual-loop open blade bent at a 90-degree angle, this retractor eliminates visual blind spots and allows surgeons to dissect, suture, and manipulate tissues directly located within the open blade aperture.

Built for stable deep-cavity clearance, this Desai Self Retaining Retractor model features an overall reach of 7.25 inches (18.41 cm, \$129). Its self-retaining chassis incorporates a multi-position mechanical ratchet lock that maintains constant, steady counter-traction on wound margins, eliminating manual assistant fatigue. The smooth wire-like perimeter of the open blades prevents focal crush trauma to delicate visceral walls and pelvic musculature during lengthy surgical procedures.

Forged from Premium Grade AISI 304/316 Stainless Steel (Non-Magnetic) with a reflection-free satin matte finish, Retractor Maker manufactures the Desai Self Retaining Retractor to withstand repetitive hospital sterilization cycles, providing general, urological, and gynecological surgical teams with high corrosion resistance, unyielding structural stability, and complete autoclave durability.



THANK YOU

We look forward to discussing your surgical retractor requirements.

CONTACT RETRACTOR MAKER

EMAIL

info@retractormaker.com

WHATSAPP

+92 327 8345101

WEBSITES

retractormaker.com
retractormanufacturer.com

ADDRESS

Retractor Maker, Near FFRV+F29, Roras Road, Sialkot-51310, Pakistan

FOLLOW US



[LinkedIn](#)



[Facebook](#)



[Instagram](#)



[X \(Twitter\)](#)



[YouTube](#)



[Pinterest](#)