



RETRACTOR MAKER

PRODUCT BROCHURE OF

Kelly Sims Vaginal Retractor

CERTIFIED & APPROVED

✓ EN ISO 13485:2016 + A11:2021



ISO 9001:2015



EN ISO 7153-1:2016



cGMP Approved




SCCI Registered since 2023



SIMAP Approved since 2023

All procedures and processes follow FDA guidelines.



Kelly Sims Vaginal Retractor

SKU RM-KSVR-1

\$84 USD

AVAILABLE SIZES & PRICING

SIZE / OPTION	SKU / DETAILS	PRICE (USD)
18 cm	RM-KSVR-1-1	\$84
20 cm	RM-KSVR-1-2	\$89
23 cm	RM-KSVR-1-3	\$94
26 cm	RM-KSVR-1-4	\$99

PRODUCT DESCRIPTION

The Kelly Sims Vaginal Retractor is an essential gynecological and pelvic reconstructive surgical instrument engineered specifically to provide deep, stable vaginal wall displacement and apical vault visualization. Handcrafted from Premium Grade AISI 304/316 Stainless Steel (Non-Magnetic), this retractor allows surgeons to depress the posterior vaginal wall and pelvic floor smoothly, establishing an expansive visual corridor to the cervix and fornices without causing focal mucosal tears.

Designed to accommodate variable anatomical depths, the Kelly Sims Vaginal Retractor can be customized across four overall lengths: 18 cm (\$84), 20 cm (\$89), 23 cm (\$94), and 26 cm (\$99). Surgeons can pair these lengths with three distinct blade footprints: 57 mm × 25 mm, 70 mm × 32 mm, and an expansive 76 mm × 38 mm surface area. The smoothly contoured, rounded blade lip prevents tissue trauma, while the ergonomic drop handle provides operating assistants with steady, non-slip manual leverage during long pelvic floor repairs.

Forged with an anti-glare satin matte finish, Retractor Maker manufactures the Kelly Sims Vaginal Retractor to withstand repetitive hospital sterilization cycles, providing obstetrical, gynecological, and urogynecological teams worldwide with exceptional durability, corrosion resistance, and complete autoclave compatibility.



THANK YOU

We look forward to discussing your surgical retractor requirements.

CONTACT RETRACTOR MAKER

EMAIL

info@retractormaker.com

WHATSAPP

+92 327 8345101

WEBSITES

retractormaker.com
retractormanufacturer.com

ADDRESS

Retractor Maker, Near FFRV+F29, Roras Road, Sialkot-51310, Pakistan

FOLLOW US



[LinkedIn](#)



[Facebook](#)



[Instagram](#)



[X \(Twitter\)](#)



[YouTube](#)



[Pinterest](#)