



RETRACTOR MAKER

PRODUCT BROCHURE OF

Ferreira breast augmentation Retractor

CERTIFIED & APPROVED

✓ EN ISO 13485:2016 + A11:2021



ISO 9001:2015



EN ISO 7153-1:2016



cGMP Approved




SCCI Registered since 2023



SIMAP Approved since 2023

All procedures and processes follow FDA guidelines.



Ferreira breast augmentation Retractor

SKU RM-FBAR-1

\$99 USD

AVAILABLE SIZES & PRICING

SIZE / OPTION	SKU / DETAILS	PRICE (USD)
Ferreira breast augmentation Retractor	5 ⅞" (130 mm) handle, 25 × 171 mm curved tapered blade	\$99
Protective Sleeve / Box	Individual transit packaging	\$99

PRODUCT DESCRIPTION

The Ferreira breast augmentation Retractor is an internationally acclaimed aesthetic surgery instrument engineered specifically to provide broad, atraumatic tissue elevation and pocket visualization during cosmetic and reconstructive breast procedures. Designed with an elongated curved blade profile, this retractor allows surgeons to dissect and elevate the pectoralis major muscle smoothly, establishing an expansive working pocket for silicone or saline implant insertion.

Built for optimal tissue clearance and safety, this Ferreira breast augmentation Retractor configuration features a 5 ⅞-inch (130 mm) handle (\$99) paired with a 25 mm wide by 171 mm long blade. Its blade geometry incorporates a curved, tapered layout engineered to safely support delicate soft tissue contours without causing localized pressure necrosis or cutaneous dimpling. The ergonomic handle provides surgical assistants with steady, fatigue-free manual leverage throughout deep pocket preparation.

Forged from Premium Grade AISI 304/316 Stainless Steel (Non-Magnetic) with a reflection-free satin matte finish, Retractor Maker manufactures the Ferreira breast augmentation Retractor to withstand repetitive hospital sterilization cycles, providing plastic, aesthetic, and oncologic breast surgery teams worldwide with superior corrosion resistance, structural rigidity, and complete autoclave durability.



THANK YOU

We look forward to discussing your surgical retractor requirements.

CONTACT RETRACTOR MAKER

EMAIL

info@retractormaker.com

WHATSAPP

+92 327 8345101

WEBSITES

retractormaker.com
retractormanufacturer.com

ADDRESS

Retractor Maker, Near FFRV+F29, Roras Road, Sialkot-51310, Pakistan

FOLLOW US



[LinkedIn](#)



[Facebook](#)



[Instagram](#)



[X \(Twitter\)](#)



[YouTube](#)



[Pinterest](#)