



RETRACTOR MAKER

PRODUCT BROCHURE OF

Biggs Face Lift (Rhytidectomy) Retractor

CERTIFIED & APPROVED

✓ EN ISO 13485:2016 + A11:2021



ISO 9001:2015



EN ISO 7153-1:2016



cGMP Approved



SCCI Registered since 2023



SIMAP Approved since 2023

All procedures and processes follow FDA guidelines.



Biggs Face Lift (Rhytidectomy) Retractor

SKU RM-BFLR-1

\$79 USD

AVAILABLE SIZES & PRICING

SIZE / OPTION	SKU / DETAILS	PRICE (USD)
Biggs Face Lift (Rhytidectomy) Retractor	8" (20.3 cm) length, 5.2 cm fenestrated blade	\$79
Protective Sleeve / Box	Individual transit packaging	\$79

PRODUCT DESCRIPTION

The Biggs Face Lift (Rhytidectomy) Retractor is an internationally acclaimed aesthetic surgery instrument engineered to provide wide, atraumatic soft-tissue retraction during facial rejuvenation procedures. Designed with an expansive 5.2 cm wide blade, this instrument distributes retractive force evenly across cutaneous and SMAS flaps, preventing tissue dimpling, ischemia, or edge tearing during extensive undermining.

Built for superior operative ergonomics, this Biggs Face Lift (Rhytidectomy) Retractor model features an overall length of 8 inches (20.3 cm). Its blade integrates a specialized fenestrated architecture that reduces overall weight while allowing operating surgeons to visualize underlying tissue planes through the blade body. A gently contoured distal lip hooks securely under the skin flap, providing steady counter-traction and minimizing assistant hand fatigue during prolonged dissections.

Forged from Premium Grade AISI 304/316 Stainless Steel (Non-Magnetic) with an anti-reflective satin matte finish, Retractor Maker manufactures the Biggs Face Lift (Rhytidectomy) Retractor to strict clinical tolerances, ensuring high corrosion resistance, comfortable handling, and complete autoclave durability for plastic and aesthetic surgery centers.



THANK YOU

We look forward to discussing your surgical retractor requirements.

CONTACT RETRACTOR MAKER

EMAIL

info@retractormaker.com

WHATSAPP

+92 327 8345101

WEBSITES

retractormaker.com
retractormanufacturer.com

ADDRESS

Retractor Maker, Near FFRV+F29, Roras Road, Sialkot-51310, Pakistan

FOLLOW US



[LinkedIn](#)



[Facebook](#)



[Instagram](#)



[X \(Twitter\)](#)



[YouTube](#)



[Pinterest](#)