



RETRACTOR MAKER

PRODUCT BROCHURE OF

Beard Distal Bicep Self-Retaining Retractor

CERTIFIED & APPROVED

✓ EN ISO 13485:2016 + A11:2021



ISO 9001:2015



EN ISO 7153-1:2016



cGMP Approved



SCCI Registered since 2023



SIMAP Approved since 2023

All procedures and processes follow FDA guidelines.



Beard Distal Bicep Self-Retaining Retractor

SKU RM-BDBSRR-1

\$299 USD

AVAILABLE SIZES & PRICING

SIZE / OPTION	SKU / DETAILS	PRICE (USD)
6.375" (16.19 cm)	RM-BDBSRR-1-1	\$299
7.5" (19.05 cm)	RM-BDBSRR-1-2	\$319

PRODUCT DESCRIPTION

The Beard Distal Bicep Self-Retaining Retractor is a specialized orthopedic surgical retractor engineered specifically to provide hands-free operative access during distal biceps tendon reattachment procedures. Designed to work within confined antecubital incisions, this instrument holds back skin and soft tissue while shielding the lateral antebrachial cutaneous nerve and radial nerve branches.

Built for refined upper extremity ergonomics, the Beard Distal Bicep Self-Retaining Retractor is manufactured in two selectable lengths: 6.375 inches (\$299) and 7.5 inches (\$319). Its 16 mm wide anatomically contoured distal blade provides focused, localized exposure over the radial tuberosity without crushing surrounding tissue. The narrow, curved handle design arches away from the surgical field, maximizing the working window for drilling, suture anchor placement, and cortical button fixation.

Forged from Premium Grade AISI 304/316 Stainless Steel (Non-Magnetic) with a reflection-free satin matte finish, Retractor Maker manufactures the Beard Distal Bicep Self-Retaining Retractor to withstand repetitive hospital sterilization cycles, providing orthopedic sports medicine and upper extremity surgeons with dependable holding stability and durability.



THANK YOU

We look forward to discussing your surgical retractor requirements.

CONTACT RETRACTOR MAKER

EMAIL

info@retractormaker.com

WHATSAPP

+92 327 8345101

WEBSITES

retractormaker.com
retractormanufacturer.com

ADDRESS

Retractor Maker, Near FFRV+F29, Roras Road, Sialkot-51310, Pakistan

FOLLOW US



[LinkedIn](#)



[Facebook](#)



[Instagram](#)



[X \(Twitter\)](#)



[YouTube](#)



[Pinterest](#)