



RETRACTOR MAKER

PRODUCT BROCHURE OF

Eastman Vaginal Retractor

CERTIFIED & APPROVED

✓ EN ISO 13485:2016 + A11:2021



ISO 9001:2015



EN ISO 7153-1:2016



cGMP Approved



SCCI Registered since 2023



SIMAP Approved since 2023

All procedures and processes follow FDA guidelines.

Donald Cameron



Eastman Vaginal Retractor

SKU RM-EVR-1

\$59 USD

AVAILABLE SIZES & PRICING

SIZE / OPTION	SKU / DETAILS	PRICE (USD)
76 mm × 38 mm — 18 cm	RM-EVR-1-1	\$59
88 mm × 38 mm — 18 cm	RM-EVR-1-2	\$62
101 mm × 38 mm — 18 cm	RM-EVR-1-3	\$64

PRODUCT DESCRIPTION

The Eastman Vaginal Retractor is an essential gynecological and pelvic reconstructive surgical instrument engineered specifically to provide broad, atraumatic displacement of the anterior or posterior vaginal wall. Designed with a wide, curved concave blade terminating in a soft retaining lip, this instrument allows surgeons to maintain clear operative visualization of the cervix, vaginal vault, and pelvic floor.

Built for clinical adaptability, this Eastman Vaginal Retractor model maintains a uniform overall instrument length of 18 cm while offering three customizable blade footprints: 76 mm × 38 mm (\$59), 88 mm × 38 mm (\$62), and an extended 101 mm × 38 mm (\$64) profile for deep pelvic access. Its broad blade architecture distributes counter-traction evenly to prevent localized mucosal tearing, while the angled handle provides surgical assistants with steady, fatigue-free manual leverage.

Forged from Premium Grade AISI 304/316 Stainless Steel (Non-Magnetic) with a reflection-free satin matte finish, Retractor Maker manufactures the Eastman Vaginal Retractor to withstand repetitive hospital sterilization cycles, providing obstetrics, gynecology, and urogynecology teams worldwide with exceptional durability and clinical dependability.



THANK YOU

We look forward to discussing your surgical retractor requirements.

CONTACT RETRACTOR MAKER

EMAIL

info@retractormaker.com

WHATSAPP

+92 327 8345101

WEBSITES

retractormaker.com
retractormanufacturer.com

ADDRESS

Retractor Maker, Near FFRV+F29, Roras Road, Sialkot-51310, Pakistan

FOLLOW US



[LinkedIn](#)



[Facebook](#)



[Instagram](#)



[X \(Twitter\)](#)



[YouTube](#)



[Pinterest](#)