



# RETRACTOR MAKER

## PRODUCT BROCHURE OF

Parker Retractor- Set Of Two

### CERTIFIED & APPROVED

✓ EN ISO 13485:2016 + A11:2021



ISO 9001:2015



EN ISO 7153-1:2016



cGMP Approved



SCCI Registered since 2023



SIMAP Approved since 2023

All procedures and processes follow FDA guidelines.



## Parker Retractor- Set Of Two

**SKU** RM-PR-SET-1

**\$44 USD**

### AVAILABLE SIZES & PRICING

| SIZE / OPTION | SKU / DETAILS | PRICE (USD) |
|---------------|---------------|-------------|
| 5 in          | RM-PR-SET-1-1 | \$44        |
| 5 1/2 in      | RM-PR-SET-1-2 | \$49        |
| 6 in          | RM-PR-SET-1-3 | \$54        |
| 7 in          | RM-PR-SET-1-4 | \$59        |
| 7 1/4 in      | RM-PR-SET-1-5 | \$64        |

### PRODUCT DESCRIPTION

The Parker Retractor- Set Of Two is an essential paired surgical instrument set engineered specifically to provide safe, atraumatic wound margin displacement and soft-tissue elevation during general, orthopedic, and plastic surgical procedures. Designed with curved C-shaped blades on opposite ends of a central handle, these retractors allow surgical teams to transition effortlessly between differing tissue depths without exchanging tools.

Built for versatile clinical utility, this Parker Retractor- Set Of Two maintains a uniform double-ended design across five selectable longitudinal sizes: 5 inches (12.7 cm, \$44), 5 ½ inches (13.97 cm, \$49), 6 inches (15.24 cm, \$54), 7 inches (17.78 cm, \$59), and 7 ¼ inches (18.42 cm, \$64). Handcrafted from Premium Grade AISI 304/316 Stainless Steel (Non-Magnetic), both working ends incorporate smooth, rounded blade lips that cradle skin edges, subcutaneous fat, and muscle planes cleanly without causing localized crush trauma. The solid flat handle provides operating assistants with steady, non-slip manual control.

Forged with an anti-glare satin matte finish, Retractor Maker manufactures the Parker Retractor- Set Of Two to withstand repetitive hospital sterilization cycles, delivering high corrosion resistance, unyielding structural rigidity, and complete autoclave durability for operating rooms, trauma centers, and outpatient surgical clinics worldwide.



# THANK YOU

We look forward to discussing your surgical retractor requirements.

## CONTACT RETRACTOR MAKER

EMAIL

**[info@retractormaker.com](mailto:info@retractormaker.com)**

WHATSAPP

**+92 327 8345101**

WEBSITES

**[retractormaker.com](http://retractormaker.com)**  
**[retractormanufacturer.com](http://retractormanufacturer.com)**

ADDRESS

Retractor Maker, Near FFRV+F29, Roras Road, Sialkot-51310, Pakistan

## FOLLOW US



[LinkedIn](#)



[Facebook](#)



[Instagram](#)



[X \(Twitter\)](#)



[YouTube](#)



[Pinterest](#)