



RETRACTOR MAKER

PRODUCT BROCHURE OF

Converse Nasal retractor, 10 cm

CERTIFIED & APPROVED

✓ EN ISO 13485:2016 + A11:2021



ISO 9001:2015



EN ISO 7153-1:2016



cGMP Approved



SCCI Registered since 2023



SIMAP Approved since 2023

All procedures and processes follow FDA guidelines.



Converse Nasal retractor, 10 cm

SKU RM-CNR-10CM-1

\$14 USD

AVAILABLE SIZES & PRICING

SIZE / OPTION	SKU / DETAILS	PRICE (USD)
10 mm × 16 mm	RM-CNR-10CM-1-1	\$14
12 mm × 30 mm	RM-CNR-10CM-1-2	\$19
12 mm × 40 mm	RM-CNR-10CM-1-3	\$24
12 mm × 50 mm	RM-CNR-10CM-1-4	\$29

PRODUCT DESCRIPTION

The Converse Nasal retractor, 10 cm is a precision-crafted rhinoplasty instrument engineered specifically to provide gentle, non-traumatic soft-tissue retraction during delicate aesthetic and reconstructive nasal procedures. Designed with a shortened 10 cm overall reach, this instrument offers high mobility and fine tactile feedback, allowing surgeons to maneuver effortlessly within tight intranasal workspaces.

Built to accommodate variable anatomical depths, this Converse Nasal retractor, 10 cm series is offered in four selectable blade sizes: 10 mm × 16 mm (\$14), 12 mm × 30 mm (\$19), 12 mm × 40 mm (\$24), and 12 mm × 50 mm (\$29). The smoothly angled blade design conforms comfortably beneath the nasal skin envelope and lateral cartilages, providing steady exposure while minimizing the risk of mucosal lacerations or cartilage fractures. Its balanced flat handle ensures fatigue-free manual control for operating assistants.

Forged from Premium Grade AISI 304/316 Stainless Steel (Non-Magnetic) with a reflection-free satin matte finish, Retractor Maker manufactures the Converse Nasal retractor, 10 cm to strict clinical tolerances, delivering exceptional durability, high corrosion resistance, and complete autoclave compatibility for plastic, ENT, and reconstructive clinics worldwide.



THANK YOU

We look forward to discussing your surgical retractor requirements.

CONTACT RETRACTOR MAKER

EMAIL

info@retractormaker.com

WHATSAPP

+92 327 8345101

WEBSITES

retractormaker.com
retractormanufacturer.com

ADDRESS

Retractor Maker, Near FFRV+F29, Roras Road, Sialkot-51310, Pakistan

FOLLOW US



[LinkedIn](#)



[Facebook](#)



[Instagram](#)



[X \(Twitter\)](#)



[YouTube](#)



[Pinterest](#)