



RETRACTOR MAKER

PRODUCT BROCHURE OF

Guttman Obstetrical Retractor

CERTIFIED & APPROVED

✓ EN ISO 13485:2016 + A11:2021



ISO 9001:2015



EN ISO 7153-1:2016



cGMP Approved



SCCI Registered since 2023



SIMAP Approved since 2023

All procedures and processes follow FDA guidelines.

Donald Cameron



Guttman Obstetrical Retractor

SKU RM-GOR-1

\$94 USD

AVAILABLE SIZES & PRICING

SIZE / OPTION	SKU / DETAILS	PRICE (USD)
7.5 cm Blade Extension — 50 mm × 70 mm base with 7.5 cm extension	RM-GOR-1-1	\$94
11.0 cm Blade Extension — 50 mm × 70 mm base with 11 cm extension	RM-GOR-1-2	\$99

PRODUCT DESCRIPTION

The Guttman Obstetrical Retractor is an essential obstetrical and gynecological surgical instrument engineered specifically to provide broad, atraumatic displacement of the urinary bladder and lower abdominal wall during cesarean deliveries and open pelvic procedures. Designed with an expansive, smoothly contoured paddle blade, this instrument shields pelvic viscera safely away from the lower uterine segment, providing an unobstructed field for infant extraction.

Built with a fixed overall length framework of 24 cm, this Guttman Obstetrical Retractor features a 50 mm × 70 mm base blade surface footprint. Surgeons can customize the instrument by selecting between two blade extension lengths: a 7.5 cm blade extension (\$94) for standard pelvic anatomy or an extended 11 cm blade extension (\$99) for deep pelvic or bariatric obstetrical interventions. Its ergonomic handle provides surgical assistants with fatigue-free manual leverage and steady counter-traction.

Forged from Premium Grade AISI 304/316 Stainless Steel (Non-Magnetic) with a reflection-free satin matte finish, Retractor Maker manufactures the Guttman Obstetrical Retractor to withstand repetitive hospital sterilization cycles, providing labor, delivery, and gynecological surgical suites worldwide with exceptional durability, corrosion resistance, and clinical dependability.



THANK YOU

We look forward to discussing your surgical retractor requirements.

CONTACT RETRACTOR MAKER

EMAIL

info@retractor maker.com

WHATSAPP

+92 327 8345101

WEBSITES

retractor maker.com
retractor manufacturer.com

ADDRESS

Retractor Maker, Near FFRV+F29, Roras Road, Sialkot-51310, Pakistan

FOLLOW US



[LinkedIn](#)



[Facebook](#)



[Instagram](#)



[X \(Twitter\)](#)



[YouTube](#)



[Pinterest](#)