



# RETRACTOR MAKER

## PRODUCT BROCHURE OF

Blount Knee Retractor

### CERTIFIED & APPROVED

✓ EN ISO 13485:2016 + A11:2021



ISO 9001:2015



EN ISO 7153-1:2016



cGMP Approved



SCCI Registered since 2023



SIMAP Approved since 2023

All procedures and processes follow FDA guidelines.



## Blount Knee Retractor

**SKU** RM-BKR-1

**\$64 USD**

### AVAILABLE SIZES & PRICING

| SIZE / OPTION                              | SKU / DETAILS | PRICE (USD) |
|--------------------------------------------|---------------|-------------|
| 6 <sup>7</sup> / <sub>8</sub> " (17.5 cm)  | RM-BKR-1-1    | \$64        |
| 7" (17.8 cm)                               | RM-BKR-1-2    | \$69        |
| 10 <sup>1</sup> / <sub>4</sub> " (26.0 cm) | RM-BKR-1-3    | \$74        |
| 10 <sup>3</sup> / <sub>8</sub> " (26.4 cm) | RM-BKR-1-4    | \$79        |

### PRODUCT DESCRIPTION

The Blount Knee Retractor is an indispensable orthopedic surgical instrument engineered specifically to provide reliable leverage and soft-tissue exposure around the femoral condyles and tibial plateau during total knee arthroplasty. Designed with an anatomical beak that hooks over bone margins, this instrument allows surgical teams to retract the patellar tendon, collateral ligaments, and capsule without tissue crushing.

Built for variable anatomical depths, the Blount Knee Retractor series offers selectable length configurations including 6 <sup>7</sup>/<sub>8</sub>" (\$64), 7" (\$69), 10 <sup>1</sup>/<sub>4</sub>" (\$74), and 10 <sup>3</sup>/<sub>8</sub>" (\$79), with standard production lines accommodating a 6 <sup>3</sup>/<sub>8</sub>" structural profile. Surgeons can configure the working end with either a Narrow Short Blade (Model variant 14093) for localized access or a Wide Blade (Model variant 14094) for broad tissue displacement. Its contoured handle provides superior counter-traction control while keeping assistant hands away from bone saws.

Forged from Premium Grade AISI 304/316 Stainless Steel (Non-Magnetic) with a reflection-free satin matte finish, Retractor Maker manufactures the Blount Knee Retractor to withstand high mechanical leverage and repetitive hospital autoclave cycles, providing orthopedic joint replacement teams with unyielding clinical durability.



# THANK YOU

We look forward to discussing your surgical retractor requirements.

## CONTACT RETRACTOR MAKER

EMAIL

**[info@retractormaker.com](mailto:info@retractormaker.com)**

WHATSAPP

**+92 327 8345101**

WEBSITES

**[retractormaker.com](http://retractormaker.com)**  
**[retractormanufacturer.com](http://retractormanufacturer.com)**

ADDRESS

Retractor Maker, Near FFRV+F29, Roras Road, Sialkot-51310, Pakistan

## FOLLOW US



[LinkedIn](#)



[Facebook](#)



[Instagram](#)



[X \(Twitter\)](#)



[YouTube](#)



[Pinterest](#)