



RETRACTOR MAKER

PRODUCT BROCHURE OF

Adson Cerebellum Retractor

CERTIFIED & APPROVED

✓ EN ISO 13485:2016 + A11:2021



ISO 9001:2015



EN ISO 7153-1:2016



cGMP Approved



SCCI Registered since 2023



SIMAP Approved since 2023

All procedures and processes follow FDA guidelines.



Adson Cerebellum Retractor

SKU RM-ACR-1

\$159 USD

AVAILABLE SIZES & PRICING

SIZE / OPTION	SKU / DETAILS	PRICE (USD)
Small — 6.5 in, 22 mm deep	RM-ACR-1-1	\$149
Standard — 7.5 in, 22 mm deep	RM-ACR-1-2	\$159
Large — 9.5 in, 22 mm deep	RM-ACR-1-3	\$174

PRODUCT DESCRIPTION

The Adson Cerebellum Retractor is an essential self-retaining surgical instrument engineered to provide steady, hands-free scalp and muscle retraction during posterior fossa cranial procedures. By distributing continuous counter-traction across heavy suboccipital musculature, this retractor eliminates assistant strain while establishing an uncompromised operative corridor for neurosurgical intervention.

Designed for anatomical ergonomics around the cranium, this Adson Cerebellum Retractor model features an overall length of 7 ½ inches (19.05 cm) with an angled arm profile that contours comfortably to patient anatomy. The working ends are equipped with sharp interlocking prongs that grip deep musculature and galea firmly, while the integrated finger-ring ratchet mechanism locks securely at variable expansion widths to maintain retraction under high muscle tension.

Manufactured from Premium Grade AISI 304/316 Stainless Steel (Non-Magnetic) with a non-reflective satin matte finish, Retractor Maker produces the Adson Cerebellum Retractor to strict clinical tolerances, ensuring smooth ratchet engagement, high corrosion resistance, and complete autoclave compatibility for neurosurgical operating rooms worldwide.



THANK YOU

We look forward to discussing your surgical retractor requirements.

CONTACT RETRACTOR MAKER

EMAIL

info@retractormaker.com

WHATSAPP

+92 327 8345101

WEBSITES

retractormaker.com
retractormanufacturer.com

ADDRESS

Retractor Maker, Near FFRV+F29, Roras Road, Sialkot-51310, Pakistan

FOLLOW US



[LinkedIn](#)



[Facebook](#)



[Instagram](#)



[X \(Twitter\)](#)



[YouTube](#)



[Pinterest](#)