



RETRACTOR MAKER

PRODUCT BROCHURE OF

Gosset Self Retaining Retractor

CERTIFIED & APPROVED

✓ EN ISO 13485:2016 + A11:2021



ISO 9001:2015



EN ISO 7153-1:2016



cGMP Approved



SCCI Registered since 2023



SIMAP Approved since 2023

All procedures and processes follow FDA guidelines.



Gosset Self Retaining Retractor

SKU RM-GSR-SET-1

\$349 USD

COMPONENT LIST

QUANTITY	COMPONENT	DIMENSIONS / DESCRIPTION
1	Self-Retaining Sliding Track Frame	Base chassis available with straight or curved mechanical shanks
2 (1 Pair)	Fixed Wire Lateral Blades	63 mm × 35 mm loop-style non-trauma blades
1	Detachable Solid Central Blade	60 mm × 65 mm slightly curved tissue retraction paddle (in 3-blade setup)
1	Dual Sliding Screw Lock Gear	Multi-point manual thumb locking screw system to anchor track positioning
2	Ergonomic Forceps Extension Shanks	Lightweight handle rails contoured with completely smooth outer edges
1	Protective Storage / Box	Standard transit packaging

PRODUCT DESCRIPTION

The Gosset Self Retaining Retractor is an internationally acclaimed abdominal surgical retractor engineered specifically to provide dependable, hands-free wound margin separation and visceral displacement during laparotomy and pelvic procedures. Designed with an adjustable slide bar chassis and non-traumatic loop-style wire lateral blades, this instrument maintains constant counter-traction across the abdominal wall without causing focal ischemic tissue contusions.

Built to clinical framework specifications, this Gosset Self Retaining Retractor complete set (\$349) features a self-retaining sliding track frame available with straight or curved mechanical shanks paired with ergonomic forceps extension shanks contoured with smooth outer edges. Handcrafted from Premium Grade AISI 304/316

Stainless Steel (Non-Magnetic), the inventory can be configured into two distinct paths: a 3-Blade Full System featuring a 60 mm × 65 mm solid central blade paired with 63 mm × 35 mm wire lateral blades expanding up to a 160 mm maximum track spread, or a 2-Blade Narrow System featuring 63 mm × 35 mm wire lateral blades expanding up to a 140 mm maximum track spread. The dual sliding screw lock gear secures position instantly.

Forged with an anti-glare satin matte finish, Retractor Maker manufactures the Gosset Self Retaining Retractor to withstand repetitive hospital autoclave cycles, providing general surgery suites, obstetrical theaters, and emergency laparotomy clinics worldwide with exceptional mechanical stability, corrosion resistance, and complete autoclave durability.



THANK YOU

We look forward to discussing your surgical retractor requirements.

CONTACT RETRACTOR MAKER

EMAIL

info@retractormaker.com

WHATSAPP

+92 327 8345101

WEBSITES

retractormaker.com
retractormanufacturer.com

ADDRESS

Retractor Maker, Near FFRV+F29, Roras Road, Sialkot-51310, Pakistan

FOLLOW US



[LinkedIn](#)



[Facebook](#)



[Instagram](#)



[X \(Twitter\)](#)



[YouTube](#)



[Pinterest](#)