



RETRACTOR MAKER

PRODUCT BROCHURE OF

APC Hip Retractor

CERTIFIED & APPROVED

✓ EN ISO 13485:2016 + A11:2021



ISO 9001:2015



EN ISO 7153-1:2016



cGMP Approved



SCCI Registered since 2023



SIMAP Approved since 2023

All procedures and processes follow FDA guidelines.



APC Hip Retractor

SKU RM-APCHR-1

\$49 USD

COMPONENT LIST

QUANTITY	COMPONENT	DIMENSIONS / DESCRIPTION
1	Single Prong	10" (25.4 cm) — targeted osseous purchase on the acetabular rim
1	Double Prongs	10" (25.4 cm) — balanced dual-point soft tissue and capsule clearance
1	Protective Sleeve / Box	Standard transit packaging

PRODUCT DESCRIPTION

The APC Hip Retractor is a heavy-duty orthopedic surgical retractor engineered specifically to provide clear, stable exposure of deep pelvic and acetabular structures during total hip replacement procedures. Designed to navigate confined joint spaces, this instrument allows surgical assistants to maintain continuous counter-traction on dense muscle and capsular layers without slipping off bony margins.

Engineered for precision leverage, this APC Hip Retractor model features an overall length of 10 inches (25.4 cm) with an angled and contoured blade profile optimized for the anatomy of the hip capsule. Surgeons can choose between a sharp single-prong tip for targeted bone anchorage or a double-prong configuration for broader soft-tissue clearance. The instrument is formatted with an extended hexagonal handle that provides superior rotational control and reduces grip fatigue during long joint reconstruction tasks.

Manufactured from Premium Grade AISI 304/316 Stainless Steel (Non-Magnetic) with a reflection-free satin matte finish, Retractor Maker produces the APC Hip Retractor to endure substantial mechanical torque and repetitive hospital autoclaving cycles, giving joint reconstruction teams reliable durability and performance.



THANK YOU

We look forward to discussing your surgical retractor requirements.

CONTACT RETRACTOR MAKER

EMAIL

info@retractormaker.com

WHATSAPP

+92 327 8345101

WEBSITES

retractormaker.com
retractormanufacturer.com

ADDRESS

Retractor Maker, Near FFRV+F29, Roras Road, Sialkot-51310, Pakistan

FOLLOW US



[LinkedIn](#)



[Facebook](#)



[Instagram](#)



[X \(Twitter\)](#)



[YouTube](#)



[Pinterest](#)