



RETRACTOR MAKER

PRODUCT BROCHURE OF

Cooley Sternotomy Retractor

CERTIFIED & APPROVED

✓ EN ISO 13485:2016 + A11:2021



ISO 9001:2015



EN ISO 7153-1:2016



cGMP Approved



SCCI Registered since 2023



SIMAP Approved since 2023

All procedures and processes follow FDA guidelines.



Cooley Sternotomy Retractor

SKU RM-CSTR-PED-1

\$359 USD

AVAILABLE SIZES & PRICING

SIZE / OPTION	SKU / DETAILS	PRICE (USD)
Newborn Profile — 85.0 mm × 55.0 mm × 16.0 mm × 14.0 mm	RM-CSTR-PED-1-1	\$359
Pediatric Profile — 115.0 mm × 80.0 mm × 30.0 mm × 15.0 mm	RM-CSTR-PED-1-2	\$369

PRODUCT DESCRIPTION

The Cooley Sternotomy Retractor is a precision-engineered pediatric cardiothoracic surgical instrument designed specifically to provide gentle, controlled sternal retraction in neonates and young children. Scaled down to accommodate smaller thoracic dimensions, this retractor establishes a wide, stable operative window into the cardiac chambers while preventing excessive focal torque on fragile infant bones.

Built for micro-mechanical accuracy, this Cooley Sternotomy Retractor series is offered in two specialized clinical profiles: a Newborn Profile (85.0 × 55.0 × 16.0 × 14.0 mm, \$359) and a Pediatric Profile (115.0 × 80.0 × 30.0 × 15.0 mm, \$369). Its finely calibrated gear rack allows millimeter-by-millimeter spread control, ensuring surgeons can achieve adequate visualization for congenital cardiac repairs without causing costochondral fractures. The contoured, polished blade edges seat smoothly against the sternal margins without slipping.

Manufactured from Premium Grade AISI 304/316 Stainless Steel (Non-Magnetic) with a reflection-free satin matte finish, Retractor Maker produces the Cooley Sternotomy Retractor to withstand repetitive hospital sterilization cycles, providing pediatric cardiac and congenital heart centers worldwide with high corrosion resistance, unyielding reliability, and clinical longevity.



THANK YOU

We look forward to discussing your surgical retractor requirements.

CONTACT RETRACTOR MAKER

EMAIL

info@retractormaker.com

WHATSAPP

+92 327 8345101

WEBSITES

retractormaker.com
retractormanufacturer.com

ADDRESS

Retractor Maker, Near FFRV+F29, Roras Road, Sialkot-51310, Pakistan

FOLLOW US



[LinkedIn](#)



[Facebook](#)



[Instagram](#)



[X \(Twitter\)](#)



[YouTube](#)



[Pinterest](#)