



RETRACTOR MAKER

PRODUCT BROCHURE OF

Denis-Browne Retractor

CERTIFIED & APPROVED

✓ EN ISO 13485:2016 + A11:2021



ISO 9001:2015



EN ISO 7153-1:2016



cGMP Approved




SCCI Registered since 2023



SIMAP Approved since 2023

All procedures and processes follow FDA guidelines.

Donald Cameron



Denis-Browne Retractor

SKU RM-DBR-SET-1

\$499 USD

COMPONENT LIST

QUANTITY	COMPONENT	DIMENSIONS / DESCRIPTION
1	Self-Retaining Circular Base Frame	175 mm × 150 mm chassis allowing all-direction visibility
2 (1 Pair)	Interchangeable Shallow Tissue Blades	40 mm × 30 mm curved tips ideal for delicate layers
2 (1 Pair)	Interchangeable Deep Tissue Blades	40 mm × 40 mm curved tips optimized for bulky tissues
4	Quick-Slide Interchanging Pin Locks	Slide-and-lock anchoring mechanism to lower finger strain
1	Protective Storage / Box	Standard transit packaging

PRODUCT DESCRIPTION

The Denis-Browne Retractor is an internationally acclaimed self-retaining surgical retractor system engineered specifically to provide expansive, 360-degree radial soft-tissue and wound margin displacement during pediatric, abdominal, and pelvic operations. Designed with a perimeter ring frame, this instrument allows surgeons to position multiple interchangeable blades at any vector along the incision, eliminating assistant hand fatigue while preserving an open, static visual corridor.

Built to clinical framework specifications, this Denis-Browne Retractor complete set (\$499) features a self-retaining circular base frame with a 175 mm × 150 mm structural ring layout that ensures all-direction visibility across deep operative fields. Handcrafted from Premium Grade AISI 304/316 Stainless Steel (Non-Magnetic), the system includes an integrated blade variety assortment containing two interchangeable shallow tissue blades (40 mm × 30 mm curved tips ideal for delicate layers) and two interchangeable deep tissue blades (40 mm × 40 mm curved tips optimized for bulky tissues). Four quick-slide interchanging pin locks

secure blades instantly to lower finger strain and speed intraoperative adjustments.

Forged with an anti-glare satin matte finish, Retractor Maker manufactures the Denis-Browne Retractor to withstand repetitive hospital autoclave cycles, providing pediatric surgical suites, urological operating theaters, and general laparotomy centers worldwide with exceptional durability, corrosion resistance, and complete autoclave durability.



THANK YOU

We look forward to discussing your surgical retractor requirements.

CONTACT RETRACTOR MAKER

EMAIL

info@retractormaker.com

WHATSAPP

+92 327 8345101

WEBSITES

retractormaker.com
retractormanufacturer.com

ADDRESS

Retractor Maker, Near FFRV+F29, Roras Road, Sialkot-51310, Pakistan

FOLLOW US



[LinkedIn](#)



[Facebook](#)



[Instagram](#)



[X \(Twitter\)](#)



[YouTube](#)



[Pinterest](#)