



RETRACTOR MAKER

PRODUCT BROCHURE OF

Cobra Retractor

CERTIFIED & APPROVED

✓ EN ISO 13485:2016 + A11:2021



ISO 9001:2015



EN ISO 7153-1:2016



cGMP Approved



SCCI Registered since 2023



SIMAP Approved since 2023

All procedures and processes follow FDA guidelines.



Cobra Retractor

SKU RM-COBR-1

\$79 USD

AVAILABLE SIZES & PRICING

SIZE / OPTION	SKU / DETAILS	PRICE (USD)
Standard Tip — 12.0"	RM-COBR-1-1	\$79
Standard Sharp Tip — 12.0"	RM-COBR-1-2	\$84
Wide Model — 11.75"	RM-COBR-1-3	\$89
Modified Small — 12.5"	RM-COBR-1-4	\$94
Double Ended — 14.0"	RM-COBR-1-5	\$99
Deep Model — 14.5"	RM-COBR-1-6	\$104
Modified Large — 14.75"	RM-COBR-1-7	\$109

PRODUCT DESCRIPTION

The Cobra Retractor is a premier orthopedic surgical instrument engineered specifically to provide powerful, stable bone leverage and deep soft-tissue retraction during total hip replacement and pelvic reconstruction surgery. Designed with an iconic curved, cobra-like flared blade that tapers to a contoured tip, this instrument hooks securely over the acetabular rim or femoral cortex, displacing large muscle masses away from the joint.

Built for complete procedural adaptability, the Cobra Retractor series is manufactured across seven specialized configurations spanning lengths from 11.75 inches to 14.75 inches: Standard Tip (12", \$79), Standard with Sharp Tip (12", \$84), Wide (11.75", \$89), Modified Small (12.5", \$94), Double Ended (14", \$99), Deep (14.5", \$104), and Modified Large (14.75", \$109). Its curved shank drops the assistant's hands safely away from the operative field, preserving direct line of sight during acetabular reaming and femoral broaching.

Forged from Premium Grade AISI 304/316 Stainless Steel (Non-Magnetic) with an anti-glare satin matte finish, Retractor Maker manufactures the Cobra Retractor to withstand immense mechanical leverage without bending, providing joint replacement and trauma surgeons worldwide with unyielding clinical durability.



THANK YOU

We look forward to discussing your surgical retractor requirements.

CONTACT RETRACTOR MAKER

EMAIL

info@retractormaker.com

WHATSAPP

+92 327 8345101

WEBSITES

retractormaker.com
retractormanufacturer.com

ADDRESS

Retractor Maker, Near FFRV+F29, Roras Road, Sialkot-51310, Pakistan

FOLLOW US



[LinkedIn](#)



[Facebook](#)



[Instagram](#)



[X \(Twitter\)](#)



[YouTube](#)



[Pinterest](#)