



RETRACTOR MAKER

PRODUCT BROCHURE OF

Terry Nasal Retractor

CERTIFIED & APPROVED

✓ EN ISO 13485:2016 + A11:2021



ISO 9001:2015



EN ISO 7153-1:2016



cGMP Approved



SCCI Registered since 2023



SIMAP Approved since 2023

All procedures and processes follow FDA guidelines.



Terry Nasal Retractor

SKU RM-TNR-6IN-1

\$44 USD

AVAILABLE SIZES & PRICING

SIZE / OPTION	SKU / DETAILS	PRICE (USD)
Terry Nasal Retractor	6" (152 mm) overall length, 5 × 46 mm blade, ergonomic handle	\$44
Protective Sleeve / Box	Individual transit packaging	\$44

PRODUCT DESCRIPTION

The Terry Nasal Retractor is a specialized handheld rhinoplasty and facial plastic surgical retractor engineered specifically to provide safe, non-traumatic elevation of the nasal skin envelope and lateral cartilages. Designed with an elongated slender blade, this instrument allows surgeons to navigate narrow subcutaneous tunnels across the nasal dorsum cleanly, shielding delicate cutaneous tissues from focal bruising or mucosal tearing during bone and cartilage remodeling.

Built to withstand high-frequency autoclave sterilization runs while offering a lifetime warranty against rust and manufacturing defects, this Terry Nasal Retractor configuration features an overall tool length of 6 inches (152 mm, \$44). Handcrafted from Premium Grade AISI 304/316 Stainless Steel (Non-Magnetic), the working blade dimensions measure 5 mm wide × 46 mm long, contoured with smooth, rounded atraumatic borders to support nasal tissue flaps safely. Its structural geometry features a single-ended blade chassis integrated into a smooth, ergonomically balanced handle that delivers fatigue-free manual leverage for operating assistants.

Provided as a non-sterile, re-sterilizable operating room device with an anti-glare satin matte finish, Retractor Maker produces the Terry Nasal Retractor to strict clinical standards, delivering high corrosion resistance, unyielding structural rigidity, and complete autoclave durability for cosmetic, ENT, and facial plastic surgery centers worldwide.



THANK YOU

We look forward to discussing your surgical retractor requirements.

CONTACT RETRACTOR MAKER

EMAIL

info@retractormaker.com

WHATSAPP

+92 327 8345101

WEBSITES

retractormaker.com
retractormanufacturer.com

ADDRESS

Retractor Maker, Near FFRV+F29, Roras Road, Sialkot-51310, Pakistan

FOLLOW US



[LinkedIn](#)



[Facebook](#)



[Instagram](#)



[X \(Twitter\)](#)



[YouTube](#)



[Pinterest](#)