



RETRACTOR MAKER

PRODUCT BROCHURE OF

Freeman Flap Retractor With Fiber Optic Illumination

CERTIFIED & APPROVED

✓ EN ISO 13485:2016 + A11:2021



ISO 9001:2015



EN ISO 7153-1:2016



cGMP Approved



SCCI Registered since 2023



SIMAP Approved since 2023

All procedures and processes follow FDA guidelines.



Freeman Flap Retractor With Fiber Optic Illumination

SKU RM-FFR-FO-1

\$89 USD

AVAILABLE SIZES & PRICING

SIZE / OPTION	SKU / DETAILS	PRICE (USD)
13 mm Wide Blade — 8" overall length, 13 mm blade, fiber-optic bundle	RM-FFR-FO-1-1	\$89
25 mm Wide Blade — 8" overall length, 25 mm blade, fiber-optic bundle	RM-FFR-FO-1-2	\$94

PRODUCT DESCRIPTION

The Freeman Flap Retractor With Fiber Optic Illumination is a premier illuminated plastic surgery retractor engineered specifically to provide superior pocket visualization and gentle flap elevation during aesthetic and reconstructive procedures. By channeling high-intensity fiber-optic light along its blade undersurface, this instrument eliminates shadowing within deep subcutaneous pockets, allowing surgeons to verify hemostasis and anatomical landmarks with ease.

Designed for precise reach and cavity clearance, this Freeman Flap Retractor With Fiber Optic Illumination model features an overall tool profile of 8 inches (203 mm) and a matching 8-inch (20.3 cm) handle. The retractor is available in two selectable blade widths: a 13 mm width (\$89) for narrow dissections and a 25 mm width (\$94) for expansive tissue flaps. The integrated optical bundle incorporates a universal fitting compatible with standard hospital light cables (Storz, Olympus, Wolf, ACMI).

Crafted from Premium Grade AISI 304/316 Stainless Steel (Non-Magnetic) and high-transmission optical glass fibers with an anti-glare satin matte finish, Retractor Maker manufactures the Freeman Flap Retractor With Fiber Optic Illumination to strict clinical standards, delivering exceptional durability, optical transmission, and complete autoclave compatibility for plastic and cosmetic surgery centers worldwide.



THANK YOU

We look forward to discussing your surgical retractor requirements.

CONTACT RETRACTOR MAKER

EMAIL

info@retractormaker.com

WHATSAPP

+92 327 8345101

WEBSITES

retractormaker.com
retractormanufacturer.com

ADDRESS

Retractor Maker, Near FFRV+F29, Roras Road, Sialkot-51310, Pakistan

FOLLOW US



[LinkedIn](#)



[Facebook](#)



[Instagram](#)



[X \(Twitter\)](#)



[YouTube](#)



[Pinterest](#)