



RETRACTOR MAKER

PRODUCT BROCHURE OF

Jackson Tracheal Retractor

CERTIFIED & APPROVED

✓ EN ISO 13485:2016 + A11:2021



ISO 9001:2015



EN ISO 7153-1:2016



cGMP Approved



SCCI Registered since 2023



SIMAP Approved since 2023

All procedures and processes follow FDA guidelines.



Jackson Tracheal Retractor

SKU RM-JTR-14-1

\$29 USD

AVAILABLE SIZES & PRICING

SIZE / OPTION	SKU / DETAILS	PRICE (USD)
Jackson Tracheal Retractor	14 cm overall length, 8.0 mm angled sharp tip, flat handle	\$29
Protective Sleeve / Box	Individual transit packaging	\$29

PRODUCT DESCRIPTION

The Jackson Tracheal Retractor is a classic airway surgical instrument engineered specifically to provide secure, non-slip elevation and stabilization of the tracheal cartilages during emergent and elective tracheostomy procedures. Designed with an angled sharp hook tip, this instrument allows surgeons to engage and elevate the cricoid or tracheal rings cleanly while displacing overlying pretracheal fascia and thyroid isthmus tissues.

Built for focused, space-limiting surgical fields, this Jackson Tracheal Retractor configuration features an overall tool length of 14 cm (\$29). Its distal working end is equipped with an 8.0 mm angled, sharp pointed tip specifically engineered to bite securely into fibrous perichondrium, preventing accidental instrument displacement during crucial airway incision and cannula insertion. The ergonomic flat-surfaced handle frame provides surgical assistants with steady, non-slip manual grip control throughout the intervention.

Forged from Premium Grade AISI 304/316 Stainless Steel (Non-Magnetic) with a reflection-free satin matte finish, Retractor Maker produces the Jackson Tracheal Retractor to withstand repetitive hospital sterilization cycles, providing ENT, trauma, and critical care operating suites worldwide with exceptional corrosion resistance, structural rigidity, and complete autoclave durability.



THANK YOU

We look forward to discussing your surgical retractor requirements.

CONTACT RETRACTOR MAKER

EMAIL

info@retractor maker.com

WHATSAPP

+92 327 8345101

WEBSITES

retractor maker.com
retractor manufacturer.com

ADDRESS

Retractor Maker, Near FFRV+F29, Roras Road, Sialkot-51310, Pakistan

FOLLOW US



[LinkedIn](#)



[Facebook](#)



[Instagram](#)



[X \(Twitter\)](#)



[YouTube](#)



[Pinterest](#)