



RETRACTOR MAKER

PRODUCT BROCHURE OF

Small Incision Accessory Set

CERTIFIED & APPROVED

✓ EN ISO 13485:2016 + A11:2021



ISO 9001:2015



EN ISO 7153-1:2016



cGMP Approved




SCCI Registered since 2023



SIMAP Approved since 2023

All procedures and processes follow FDA guidelines.



Small Incision Accessory Set

SKU RM-SIAS-1

\$449 USD

COMPONENT LIST

QUANTITY	COMPONENT	DIMENSIONS / DESCRIPTION
1	Gelpi Hook Blade	Standard Fine-Point Profile
1	Kelly Blade	1" × 1½" (2.5 × 3.8 cm)
1	Kelly Blade	1" × 2½" (2.5 × 6.3 cm)
1	Small Round Solid Ring	8.5" (21.3 cm) Inner Diameter (ID)
1	Sterilization Box	Dedicated Storage Organizer

PRODUCT DESCRIPTION

The Small Incision Accessory Set is an add-on kit for surgical teams that already own a Bookwalter General Surgery Retractor System and work through limited exposures. It is not a standalone system — it is fully compatible with Bookwalter General Surgery Retractor Systems and standard Bookwalter-style components.

The kit contains a Gelpi hook blade with a standard fine-point profile for precise, narrow retraction, along with 1" × 1½" (2.5 × 3.8 cm) and 1" × 2½" (2.5 × 6.3 cm) Kelly blades for shallow retraction at two depths. An 8.5" (21.3 cm) inner-diameter small round solid ring gives a compact frame that keeps hardware clear of a small incision.

All components are manufactured from Premium Grade AISI 304/316 Stainless Steel (Non-Magnetic) and supplied in a dedicated sterilization box. The Small Incision Accessory Set is supplied non-sterile, is reusable and fully autoclavable, and must be cleaned and sterilized according to your facility's validated procedures before use.

Retractor Maker manufactures and supplies the set directly from Sialkot, Pakistan. The Small Incision Accessory Set is manufactured, inspected and supplied directly by Retractor Maker.



THANK YOU

We look forward to discussing your surgical retractor requirements.

CONTACT RETRACTOR MAKER

EMAIL

info@retractormaker.com

WHATSAPP

+92 327 8345101

WEBSITES

retractormaker.com
retractormanufacturer.com

ADDRESS

Retractor Maker, Near FFRV+F29, Roras Road, Sialkot-51310, Pakistan

FOLLOW US



[LinkedIn](#)



[Facebook](#)



[Instagram](#)



[X \(Twitter\)](#)



[YouTube](#)



[Pinterest](#)