



RETRACTOR MAKER

PRODUCT BROCHURE OF

Freeman Face Lift Retractors

CERTIFIED & APPROVED

✓ EN ISO 13485:2016 + A11:2021



ISO 9001:2015



EN ISO 7153-1:2016



cGMP Approved



SCCI Registered since 2023



SIMAP Approved since 2023

All procedures and processes follow FDA guidelines.



Freeman Face Lift Retractors

SKU RM-FFLR-SERIES-1

\$62 USD

AVAILABLE SIZES & PRICING

SIZE / OPTION	SKU / DETAILS	PRICE (USD)
20 mm - 22 mm — Straight	RM-FFLR-SERIES-1-1	\$62
20 mm - 22 mm — Fully Curved	RM-FFLR-SERIES-1-2	\$64
35 mm - 36 mm — Straight	RM-FFLR-SERIES-1-3	\$66
35 mm - 36 mm — Fully Curved	RM-FFLR-SERIES-1-4	\$68
37 mm - 38 mm — Straight	RM-FFLR-SERIES-1-5	\$70
37 mm - 38 mm — Fully Curved	RM-FFLR-SERIES-1-6	\$72

PRODUCT DESCRIPTION

Freeman Face Lift Retractors are foundational aesthetic plastic surgery instruments engineered specifically to provide wide, non-slip tissue elevation and optimal visualization during rhytidectomy and facial reconstructive surgery. Designed with an expansive, contoured working end, these retractors distribute manual counter-traction evenly across skin flaps, preventing localized ischemia, dimpling, or dermal tearing during extensive undermining.

Built to suit specific anatomical zones and procedural techniques, this Freeman Face Lift Retractors line maintains a fixed overall tool length of 7 inches (17.8 cm). The series is offered across six standardized width diameters: 20 mm (\$62), 22 mm (\$64), 35 mm (\$66), 36 mm (\$68), 37 mm (\$70), and 38 mm (\$72). Surgeons can choose between a Straight Profile for direct planar lifting or a Fully Curved Pattern that arches around facial contours. Its ergonomic flat handle provides surgical assistants with steady, fatigue-free control.

Machined from Premium Grade AISI 304/316 Stainless Steel (Non-Magnetic) with a reflection-free satin matte finish, Retractor Maker manufactures Freeman Face Lift Retractors to strict clinical tolerances, delivering exceptional corrosion resistance, unyielding structural rigidity, and complete autoclave durability for cosmetic and facial plastic surgery clinics worldwide.



THANK YOU

We look forward to discussing your surgical retractor requirements.

CONTACT RETRACTOR MAKER

EMAIL

info@retractormaker.com

WHATSAPP

+92 327 8345101

WEBSITES

retractormaker.com
retractormanufacturer.com

ADDRESS

Retractor Maker, Near FFRV+F29, Roras Road, Sialkot-51310, Pakistan

FOLLOW US



[LinkedIn](#)



[Facebook](#)



[Instagram](#)



[X \(Twitter\)](#)



[YouTube](#)



[Pinterest](#)