



RETRACTOR MAKER

PRODUCT BROCHURE OF

Sheen Nasal Retractor Angled

CERTIFIED & APPROVED

✓ EN ISO 13485:2016 + A11:2021



ISO 9001:2015



EN ISO 7153-1:2016



cGMP Approved




SCCI Registered since 2023



SIMAP Approved since 2023

All procedures and processes follow FDA guidelines.



Sheen Nasal Retractor Angled

SKU RM-SNRA-65-1

\$49 USD

AVAILABLE SIZES & PRICING

SIZE / OPTION	SKU / DETAILS	PRICE (USD)
Sheen Nasal Retractor Angled	6 ½" (165 mm) overall length, 10 mm solid blade, angled chassis	\$49
Protective Sleeve / Box	Individual transit packaging	\$49

PRODUCT DESCRIPTION

The Sheen Nasal Retractor Angled is a specialized handheld rhinoplasty and facial plastic surgery instrument engineered specifically to provide safe, non-traumatic elevation of the nasal skin envelope and lower lateral cartilages. Designed with an anatomical angled shank, this instrument drops the assistant's hand safely away from the surgical axis, allowing the operating surgeon to maintain an unobstructed line of sight across the nasal dorsum.

Built to withstand frequent autoclave sterilization cycles while providing a lifetime warranty against rust, this Sheen Nasal Retractor Angled configuration features an overall tool length of 6 ½ inches (165 mm, \$49). Handcrafted from Premium Grade AISI 304/316 Stainless Steel (Non-Magnetic), the working blade width measures 10 mm wide in a solid blade layout contoured with rounded, atraumatic margins to prevent mucosal tears or dermal ischemia. Its fixed single-ended angled chassis provides operating assistants with steady, fatigue-free manual leverage throughout intricate nasal tip and vault reshaping.

Classified as a certified premium operating room grade retractor, Retractor Maker manufactures the Sheen Nasal Retractor Angled with an anti-glare satin matte finish, delivering high corrosion resistance, unyielding structural rigidity, and complete autoclave durability for cosmetic, ENT, and facial plastic surgery centers worldwide.



THANK YOU

We look forward to discussing your surgical retractor requirements.

CONTACT RETRACTOR MAKER

EMAIL

info@retractor maker.com

WHATSAPP

+92 327 8345101

WEBSITES

retractor maker.com
retractor manufacturer.com

ADDRESS

Retractor Maker, Near FFRV+F29, Roras Road, Sialkot-51310, Pakistan

FOLLOW US



[LinkedIn](#)



[Facebook](#)



[Instagram](#)



[X \(Twitter\)](#)



[YouTube](#)



[Pinterest](#)