



RETRACTOR MAKER

PRODUCT BROCHURE OF

Leaflet Retractor

CERTIFIED & APPROVED

✓ EN ISO 13485:2016 + A11:2021



ISO 9001:2015



EN ISO 7153-1:2016



cGMP Approved




SCCI Registered since 2023



SIMAP Approved since 2023

All procedures and processes follow FDA guidelines.



Leaflet Retractor

SKU RM-LFR-22-1

\$49 USD

COMPONENT LIST

QUANTITY	COMPONENT	DIMENSIONS / DESCRIPTION
1	5.0 mm Wide Blade	Overall Length: 22 cm · Primary Clinical Role: Pediatric valve surgery, commissurotomy, and aortic cusp inspection · Price (USD): \$49.00
1	8.0 mm Wide Blade	Overall Length: 22 cm · Primary Clinical Role: Adult mitral leaflet retraction and subvalvular chordal exposure · Price (USD): \$49.00
1	Protective Sleeve / Box	Overall Length: Standard Transit Packaging · Primary Clinical Role: — · Price (USD): Included

PRODUCT DESCRIPTION

The Leaflet Retractor is an indispensable cardiovascular surgical instrument engineered specifically to provide safe, non-traumatic retraction and manipulation of aortic, mitral, and tricuspid valve leaflets during open heart reconstructions. Designed with an ultra-smooth distal paddle, this instrument allows cardiothoracic surgeons to inspect subvalvular apparatus, commissures, and chordae tendineae without causing leaflet contusions or tears.

Built for micro-precision within deep cardiac chambers, this Leaflet Retractor maintains a uniform longitudinal profile of 22 cm (\$49) that keeps assistant hands well outside the intracardiac surgical field. Handcrafted from Premium Grade AISI 304/316 Stainless Steel (Non-Magnetic), the instrument is supplied with two distinct blade width options: a 5.0 mm wide blade for pediatric valves and narrow commissural corridors, and an 8.0 mm wide blade for adult aortic and mitral leaflet deflection. The balanced flat handle ensures fatigue-free manual leverage and high tactile sensitivity during intricate suturing.

Forged with an anti-glare satin matte finish, Retractor Maker manufactures the Leaflet Retractor to withstand repetitive hospital autoclave cycles, providing cardiac surgical suites worldwide with exceptional corrosion

resistance, unyielding structural rigidity, and complete clinical durability.



THANK YOU

We look forward to discussing your surgical retractor requirements.

CONTACT RETRACTOR MAKER

EMAIL

info@retractormaker.com

WHATSAPP

+92 327 8345101

WEBSITES

retractormaker.com
retractormanufacturer.com

ADDRESS

Retractor Maker, Near FFRV+F29, Roras Road, Sialkot-51310, Pakistan

FOLLOW US



[LinkedIn](#)



[Facebook](#)



[Instagram](#)



[X \(Twitter\)](#)



[YouTube](#)



[Pinterest](#)