



RETRACTOR MAKER

PRODUCT BROCHURE OF

Expanded Vascular Lobster Retractor

CERTIFIED & APPROVED

✓ EN ISO 13485:2016 + A11:2021



ISO 9001:2015



EN ISO 7153-1:2016



cGMP Approved



SCCI Registered since 2023



SIMAP Approved since 2023

All procedures and processes follow FDA guidelines.



Expanded Vascular Lobster Retractor

SKU RM-EVLR-1

\$2,999 USD

COMPONENT LIST

QUANTITY	COMPONENT	DIMENSIONS / DESCRIPTION
6	Lobster Retractor Clamp	—
1	Lobster Table Clamp with Vertical Post & Frame Coupler	—
1	Lobster Hinged Extension Arm	—
1	Lobster Frame	—
2	Lobster Mayo Swivel Retractor	2¾" × 2"
1	Lobster Mayo Swivel Retractor	3½" × 2"
1	Lobster Deaver Swivel Retractor	2¾" × 4½"
1	Lobster Harrington Sweetheart Retractor	2" × 5"
1	Lobster Malleable Swivel Retractor (Perforated)	2" × 6½"
1	Lobster Malleable Swivel Retractor (Perforated)	3" × 6½"
1	Lobster Fence Swivel Retractor	4" × 6"
2	Lobster Renal Swivel Retractor	1" × 6"
1	Lobster Splanchnic Swivel Retractor	2" × 4½"

QUANTITY	COMPONENT	DIMENSIONS / DESCRIPTION
1	Lobster Splanchnic Swivel Retractor	2" × 6"
1	Sterilization Box	—

PRODUCT DESCRIPTION

The Expanded Vascular Lobster Retractor is an extended-reach, table-mounted self-retaining retractor engineered to meet the stringent visualization requirements of complex vascular reconstructions. By combining an articulated hinged extension arm with a rigid Lobster frame and table post coupler, this system provides unmatched reach into deep retroperitoneal fields while eliminating surgical assistant fatigue.

Designed for complicated aortic and peripheral vessel exposure, the Expanded Vascular Lobster Retractor configuration incorporates six quick-locking slide clamps and an expanded 12-piece swivel blade set. The complement includes dual narrow Renal blades, multiple deep visceral Splanchnic retractors, broad protective Fence blades, Harrington Sweethearts, Deaver, Mayo, and bendable perforated malleable ribbons that rotate 360 degrees to conform atraumatically to vascular margins.

Machined from Premium Grade AISI 304/316 Stainless Steel (Non-Magnetic) and high-strength alloys with an anti-glare satin finish, Retractor Maker supplies the complete Expanded Vascular Lobster Retractor non-sterile and fully autoclavable in an organized heavy-duty sterilization box, providing surgical teams with maximum clinical versatility.



THANK YOU

We look forward to discussing your surgical retractor requirements.

CONTACT RETRACTOR MAKER

EMAIL

info@retractormaker.com

WHATSAPP

+92 327 8345101

WEBSITES

retractormaker.com
retractormanufacturer.com

ADDRESS

Retractor Maker, Near FFRV+F29, Roras Road, Sialkot-51310, Pakistan

FOLLOW US



[LinkedIn](#)



[Facebook](#)



[Instagram](#)



[X \(Twitter\)](#)



[YouTube](#)



[Pinterest](#)