



RETRACTOR MAKER

PRODUCT BROCHURE OF

Jackson Tracheal Retractor-Double Ended

CERTIFIED & APPROVED

✓ EN ISO 13485:2016 + A11:2021



ISO 9001:2015



EN ISO 7153-1:2016



cGMP Approved



SCCI Registered since 2023



SIMAP Approved since 2023

All procedures and processes follow FDA guidelines.

Donald Cameron



Jackson Tracheal Retractor-Double Ended

SKU RM-JTR-DE-17-1

\$39 USD

AVAILABLE SIZES & PRICING

SIZE / OPTION	SKU / DETAILS	PRICE (USD)
Jackson Tracheal Retractor-Double Ended	17 cm overall length, slender curved dual-ended blades	\$39
Protective Sleeve / Box	Individual transit packaging	\$39

PRODUCT DESCRIPTION

The Jackson Tracheal Retractor-Double Ended is an essential otolaryngology and emergency surgical instrument engineered specifically to provide dual-retraction efficiency during elective and emergent tracheostomy procedures. Designed with contrasting curved blades on opposite ends of a central handle, this instrument allows surgical teams to transition effortlessly between pretracheal soft-tissue displacement and direct tracheal ring elevation without exchanging tools.

Built to optimize operative agility, this Jackson Tracheal Retractor-Double Ended configuration features an overall tool length of 17 cm (\$39). Its slender, curved, double-ended architecture is paired with an ergonomic center handle designed to reduce user hand fatigue while providing steady manual counter-traction. The smoothly polished working ends cradle delicate cervical tissues and tracheal margins comfortably, preventing accidental mucosal laceration or cartilage fractures during tube insertion.

Machined from Premium Grade AISI 304/316 Stainless Steel (Non-Magnetic) with a reflection-free satin matte finish, Retractor Maker manufactures the Jackson Tracheal Retractor-Double Ended to strict clinical standards, delivering high corrosion resistance, unyielding structural rigidity, and complete autoclave durability for ENT suites, emergency rooms, and trauma centers worldwide.



THANK YOU

We look forward to discussing your surgical retractor requirements.

CONTACT RETRACTOR MAKER

EMAIL

info@retractormaker.com

WHATSAPP

+92 327 8345101

WEBSITES

retractormaker.com
retractormanufacturer.com

ADDRESS

Retractor Maker, Near FFRV+F29, Roras Road, Sialkot-51310, Pakistan

FOLLOW US



[LinkedIn](#)



[Facebook](#)



[Instagram](#)



[X \(Twitter\)](#)



[YouTube](#)



[Pinterest](#)