



RETRACTOR MAKER

PRODUCT BROCHURE OF

Rake Finger Retractor

CERTIFIED & APPROVED

✓ EN ISO 13485:2016 + A11:2021



ISO 9001:2015



EN ISO 7153-1:2016



cGMP Approved



SCCI Registered since 2023



SIMAP Approved since 2023

All procedures and processes follow FDA guidelines.



Rake Finger Retractor

SKU RM-RFR-153-1

\$29 USD

AVAILABLE SIZES & PRICING

SIZE / OPTION	SKU / DETAILS	PRICE (USD)
1 Prong (Single Hook) — Blunt or Sharp tips	RM-RFR-153-1-1	\$29
2 Prongs — Blunt or Sharp tips	RM-RFR-153-1-2	\$34
3 Prongs — Blunt or Sharp tips	RM-RFR-153-1-3	\$39
4 Prongs — Blunt or Sharp tips	RM-RFR-153-1-4	\$44

PRODUCT DESCRIPTION

The Rake Finger Retractor is a specialized handheld surgical rake retractor engineered specifically to provide precise, non-slip skin flap elevation and soft-tissue retraction during hand, plastic, and minor general surgical procedures. Designed with multiple curved rake prongs, this instrument allows surgeons to grasp cutaneous and subcutaneous margins securely without slipping or crushing delicate tissue borders.

Built for delicate manual control, this Rake Finger Retractor holds a uniform total length frame of 15.3 cm that provides operating assistants with exceptional tactile feedback and maneuverability in restricted operative fields. Handcrafted from Premium Grade AISI 304/316 Stainless Steel (Non-Magnetic), the instrument is highly customizable across four prong quantity configurations: 1 prong (\$29), 2 prongs (\$34), 3 prongs (\$39), or 4 prongs rake styles (\$44). Surgeons can select from two prong tip profiles: Blunt tips engineered for safe soft tissue manipulation minimizing trauma, or Sharp tips engineered to provide a firm, slippage-free grip on tougher tissue layers. The ergonomic flat handle ensures fatigue-free manual leverage throughout intricate dissections.

Forged with an anti-glare satin matte finish, Retractor Maker manufactures the Rake Finger Retractor to withstand repetitive hospital autoclave cycles, providing plastic, hand, and general surgery clinics worldwide

with exceptional corrosion resistance, structural rigidity, and complete clinical durability.



THANK YOU

We look forward to discussing your surgical retractor requirements.

CONTACT RETRACTOR MAKER

EMAIL

info@retractormaker.com

WHATSAPP

+92 327 8345101

WEBSITES

retractormaker.com
retractormanufacturer.com

ADDRESS

Retractor Maker, Near FFRV+F29, Roras Road, Sialkot-51310, Pakistan

FOLLOW US



[LinkedIn](#)



[Facebook](#)



[Instagram](#)



[X \(Twitter\)](#)



[YouTube](#)



[Pinterest](#)