



RETRACTOR MAKER

PRODUCT BROCHURE OF

Anterolateral Acetabular Rim Retractor

CERTIFIED & APPROVED

✓ EN ISO 13485:2016 + A11:2021



ISO 9001:2015



EN ISO 7153-1:2016



cGMP Approved



SCCI Registered since 2023



SIMAP Approved since 2023

All procedures and processes follow FDA guidelines.

Donald Cameron



Anterolateral Acetabular Rim Retractor

SKU RM-AARR-1

\$49 USD

AVAILABLE SIZES & PRICING

SIZE / OPTION	SKU / DETAILS	PRICE (USD)
Anterolateral Acetabular Rim Retractor	10" (25.4 cm) length, 18 mm width x 3.25" depth blade	\$49
Protective Sleeve / Box	Individual transit packaging	\$49

PRODUCT DESCRIPTION

The Anterolateral Acetabular Rim Retractor is a dedicated orthopedic surgical retractor engineered to establish clear, stable exposure of the acetabular margin during primary and revision total hip arthroplasty. By anchoring securely over the anterolateral bone rim, this instrument displaces anterior capsular tissues and abductor musculature, ensuring a clear visual corridor into the pelvic socket.

Designed for robust anatomical leverage, this Anterolateral Acetabular Rim Retractor model features an overall length of 10 inches (25.4 cm) with an 18 mm blade width and an extended 3.25-inch (8.3 cm) blade depth profile. The angled distal blade is contoured to hook firmly over the acetabular rim without impinging on surrounding femoral neurovascular structures. Its ergonomic handle allows surgical assistants to apply sustained counter-traction with minimal hand fatigue during socket reaming and cup impaction.

Manufactured from Premium Grade AISI 304/316 Stainless Steel (Non-Magnetic) with a glare-resistant satin matte finish, Retractor Maker produces the Anterolateral Acetabular Rim Retractor to endure heavy operative leverage, providing orthopedic departments worldwide with superior corrosion resistance and complete autoclave durability.



THANK YOU

We look forward to discussing your surgical retractor requirements.

CONTACT RETRACTOR MAKER

EMAIL

info@retractormaker.com

WHATSAPP

+92 327 8345101

WEBSITES

retractormaker.com
retractormanufacturer.com

ADDRESS

Retractor Maker, Near FFRV+F29, Roras Road, Sialkot-51310, Pakistan

FOLLOW US



[LinkedIn](#)



[Facebook](#)



[Instagram](#)



[X \(Twitter\)](#)



[YouTube](#)



[Pinterest](#)