



RETRACTOR MAKER

PRODUCT BROCHURE OF

Tebbetts Breast Retractor

CERTIFIED & APPROVED

✓ EN ISO 13485:2016 + A11:2021



ISO 9001:2015



EN ISO 7153-1:2016



cGMP Approved



SCCI Registered since 2023



SIMAP Approved since 2023

All procedures and processes follow FDA guidelines.



Tebbetts Breast Retractor

SKU RM-TBR-64-1

\$64 USD

AVAILABLE SIZES & PRICING

SIZE / OPTION	SKU / DETAILS	PRICE (USD)
9 cm × 24 mm — With Teeth or Without Teeth — Inframammary access, periareolar pocket dissection	RM-TBR-64-1-1	\$64
15 cm × 24 mm — With Teeth or Without Teeth — Deep subpectoral pocket creation and transaxillary breast augmentation	RM-TBR-64-1-2	\$64

PRODUCT DESCRIPTION

The Tebbetts Breast Retractor is an internationally acclaimed aesthetic plastic surgery retractor engineered specifically to provide precise tissue elevation and pocket visualization during breast augmentation and mammoplasty. Designed with a curved, tapered blade profile, this instrument allows surgeons to dissect and elevate the pectoralis major muscle cleanly, establishing an unhindered surgical pocket for breast prosthesis insertion.

Built to offer excellent structural durability and a lifetime warranty against rust, this Tebbetts Breast Retractor series (\$64) is highly customizable through multiple blade sizes and terminal edge configurations. Handcrafted from Premium Grade AISI 304/316 Stainless Steel (Non-Magnetic), the instrument offers two blade dimensions: a 9 cm × 24 mm layout for standard inframammary incisions, and an extended 15 cm × 24 mm layout for deep subpectoral pocket dissection. Surgeons can select between a With Teeth profile layout engineered to secure an optimal tissue grip and prevent intraoperative slippage, or a Without Teeth (Smooth End) profile layout engineered for safe, extra-traumatic soft tissue manipulation. Its ergonomic handle delivers fatigue-free manual leverage.

Forged with an anti-glare satin matte finish, Retractor Maker manufactures the Tebbetts Breast Retractor to withstand repetitive hospital autoclave cycles, providing plastic, cosmetic, and oncologic breast reconstructive suites worldwide with unyielding structural reliability, corrosion resistance, and complete clinical durability.



THANK YOU

We look forward to discussing your surgical retractor requirements.

CONTACT RETRACTOR MAKER

EMAIL

info@retractormaker.com

WHATSAPP

+92 327 8345101

WEBSITES

retractormaker.com
retractormanufacturer.com

ADDRESS

Retractor Maker, Near FFRV+F29, Roras Road, Sialkot-51310, Pakistan

FOLLOW US



[LinkedIn](#)



[Facebook](#)



[Instagram](#)



[X \(Twitter\)](#)



[YouTube](#)



[Pinterest](#)