



RETRACTOR MAKER

PRODUCT BROCHURE OF

Bent Hohmann Retractor

CERTIFIED & APPROVED

✓ EN ISO 13485:2016 + A11:2021



ISO 9001:2015



EN ISO 7153-1:2016



cGMP Approved



SCCI Registered since 2023



SIMAP Approved since 2023

All procedures and processes follow FDA guidelines.



Bent Hohmann Retractor

SKU RM-BHR-ORTHO-1

\$69 USD

AVAILABLE SIZES & PRICING

SIZE / OPTION	SKU / DETAILS	PRICE (USD)
8"	RM-BHR-ORTHO-1-1	\$69
9 ¾"	RM-BHR-ORTHO-1-2	\$74
11"	RM-BHR-ORTHO-1-3	\$79
11 ½"	RM-BHR-ORTHO-1-4	\$84
12 ⅛"	RM-BHR-ORTHO-1-5	\$89

PRODUCT DESCRIPTION

The Bent Hohmann Retractor is an indispensable orthopedic bone leverage instrument engineered to deliver exceptional tissue displacement and bone exposure during joint reconstruction and fracture stabilization. Designed with a distinct downward bend along its shank, this instrument drops the assistant's hands and handle away from the incision, preserving an unobstructed visual line of sight into deep joint spaces.

Built for high-leverage bone and muscle elevation, the Bent Hohmann Retractor series offers an extensive selection of anatomical lengths and handle profiles. Surgeons can select from overall lengths of 8" (\$69), 9 ¾" (\$74), 11" (\$79), 11 ½" (\$84), and 12 ⅛" (\$89), paired with 7", 9 ¾", or 10" handle profiles. The narrow, contoured beak seats securely behind the bone cortex to provide powerful counter-traction without slipping into neurovascular zones.

Machined from high-tensile, Premium Grade AISI 304/316 Stainless Steel (Non-Magnetic) with a reflection-free satin matte finish, Retractor Maker manufactures the Bent Hohmann Retractor to endure heavy torque without bending, providing joint replacement and trauma teams worldwide with reliable clinical performance.



THANK YOU

We look forward to discussing your surgical retractor requirements.

CONTACT RETRACTOR MAKER

EMAIL

info@retractormaker.com

WHATSAPP

+92 327 8345101

WEBSITES

retractormaker.com
retractormanufacturer.com

ADDRESS

Retractor Maker, Near FFRV+F29, Roras Road, Sialkot-51310, Pakistan

FOLLOW US



[LinkedIn](#)



[Facebook](#)



[Instagram](#)



[X \(Twitter\)](#)



[YouTube](#)



[Pinterest](#)