



RETRACTOR MAKER

PRODUCT BROCHURE OF

Goodyear Uvula Retractor

CERTIFIED & APPROVED

✓ EN ISO 13485:2016 + A11:2021



ISO 9001:2015



EN ISO 7153-1:2016



cGMP Approved



SCCI Registered since 2023



SIMAP Approved since 2023

All procedures and processes follow FDA guidelines.



Goodyear Uvula Retractor

SKU RM-GUR-18-1

\$34 USD

AVAILABLE SIZES & PRICING

SIZE / OPTION	SKU / DETAILS	PRICE (USD)
6 1/4 in (15.9 cm)	RM-GUR-18-1-1	\$26
7 3/8 in (18.7 cm) — Standard, 15 mm blade	RM-GUR-18-1-2	\$34

PRODUCT DESCRIPTION

The Goodyear Uvula Retractor is a specialized otolaryngology surgical instrument engineered specifically to provide safe, non-traumatic retraction and deflection of the uvula and soft palate during pharyngeal and palate procedures. Designed with an anatomically curved distal blade, this instrument allows ENT surgeons to displace mucosal folds smoothly, establishing an open visual field into the posterior pharynx without causing tissue bruising or uvular edema.

Built for exceptional intraoral control, this Goodyear Uvula Retractor configuration features an overall tool length of 18 cm (\$34) that ensures adequate reach and maneuverability while keeping the surgeon's hands outside the line of sight. Its working end incorporates a 15.0 mm wide smooth blade profile specifically contoured to cradle the uvula securely without focal pinching. The balanced handle provides surgical assistants with fatigue-free manual leverage throughout adenoid and palate interventions.

Forged from Premium Grade AISI 304/316 Stainless Steel (Non-Magnetic) with an anti-glare satin matte finish, Retractor Maker manufactures the Goodyear Uvula Retractor to withstand repetitive hospital sterilization cycles, providing ENT operating suites, outpatient clinics, and plastic surgery centers with high corrosion resistance, unyielding structural rigidity, and complete autoclave durability.



THANK YOU

We look forward to discussing your surgical retractor requirements.

CONTACT RETRACTOR MAKER

EMAIL

info@retractormaker.com

WHATSAPP

+92 327 8345101

WEBSITES

retractormaker.com
retractormanufacturer.com

ADDRESS

Retractor Maker, Near FFRV+F29, Roras Road, Sialkot-51310, Pakistan

FOLLOW US



[LinkedIn](#)



[Facebook](#)



[Instagram](#)



[X \(Twitter\)](#)



[YouTube](#)



[Pinterest](#)