



RETRACTOR MAKER

PRODUCT BROCHURE OF

Biggs Mammoplasty Retractor, Fiber Optics

CERTIFIED & APPROVED

✓ EN ISO 13485:2016 + A11:2021



ISO 9001:2015



EN ISO 7153-1:2016



cGMP Approved



SCCI Registered since 2023



SIMAP Approved since 2023

All procedures and processes follow FDA guidelines.



Biggs Mammaplasty Retractor, Fiber Optics

SKU RM-BMRFO-1

\$239 USD

COMPONENT LIST

QUANTITY	COMPONENT	DIMENSIONS / DESCRIPTION
1	Biggs Mammaplasty Retractor, Fiber Optics	50 mm fenestrated blade, light guide, suction tube
1	Universal Light Cable Fitting	Connects to standard OR fiber-optic cables
1	Protective Sleeve / Box	Standard transit packaging

PRODUCT DESCRIPTION

The Biggs Mammaplasty Retractor, Fiber Optics is a premier illuminated breast surgery instrument engineered to provide superior pocket visualization and active field clearance during aesthetic and reconstructive mammaplasty. By integrating high-intensity fiber-optic illumination with a continuous smoke and fluid suction channel, this retractor eliminates dark recesses and evacuates electrosurgical plume simultaneously.

Designed for demanding surgical procedures, this Biggs Mammaplasty Retractor, Fiber Optics configuration features a 50 mm standard width blade (\$239). Its working end incorporates a specialized fenestrated architecture that delivers excellent mechanical traction and stability on glandular tissues while minimizing localized contact pressure. The integrated fiber-optic bundle features a universal fitting compatible with standard hospital light cables (Storz, Olympus, Wolf, ACMI), ensuring seamless OR integration.

Crafted from Premium Grade AISI 304/316 Stainless Steel (Non-Magnetic) and high-transmission optical glass fibers with an anti-glare satin finish, Retractor Maker produces the Biggs Mammaplasty Retractor, Fiber Optics non-sterile and fully autoclavable, delivering exceptional durability, optical clarity, and clinical utility for breast surgery centers worldwide.



THANK YOU

We look forward to discussing your surgical retractor requirements.

CONTACT RETRACTOR MAKER

EMAIL

info@retractormaker.com

WHATSAPP

+92 327 8345101

WEBSITES

retractormaker.com
retractormanufacturer.com

ADDRESS

Retractor Maker, Near FFRV+F29, Roras Road, Sialkot-51310, Pakistan

FOLLOW US



[LinkedIn](#)



[Facebook](#)



[Instagram](#)



[X \(Twitter\)](#)



[YouTube](#)



[Pinterest](#)