



RETRACTOR MAKER

PRODUCT BROCHURE OF

Morse-Favaloro Rib Spreader

CERTIFIED & APPROVED

✓ EN ISO 13485:2016 + A11:2021



ISO 9001:2015



EN ISO 7153-1:2016



cGMP Approved



SCCI Registered since 2023



SIMAP Approved since 2023

All procedures and processes follow FDA guidelines.



Morse-Favaloro Rib Spreader

SKU RM-MFRS-160-1

\$269 USD

COMPONENT LIST

QUANTITY	COMPONENT	DIMENSIONS / DESCRIPTION
1	Parallel Slide Guide Rail Track Frame	Lightweight aluminum guide bar rail serving as the expansion chassis
4 (2 Pairs)	Movable Angled Lateral Blades	17 mm × 22 mm four multi-position adjustable blades that contort to ribs
1	Turning Thread Screw Lock	Integrated locking bolt that secures blades at optimum retraction
2	Ergonomic Retracting Forceps Arms	Tapered aluminum handle shanks engineered with smooth outer edges
1	Protective Storage / Box	Standard transit packaging

PRODUCT DESCRIPTION

The Morse-Favaloro Rib Spreader is an essential cardiothoracic surgical retractor engineered specifically to provide safe, contoured intercostal rib separation and anterior thoracic exposure during coronary revascularization and pulmonary procedures. Designed with multi-position movable blades, this instrument allows the lateral paddles to swivel and contort precisely to the patient's rib curvature, distributing distraction loads evenly to prevent focal periosteal bruising or costal fracturing.

Built for refined surgical control, this dedicated Morse-Favaloro Rib Spreader configuration features an overall horizontal expansion spread of 160 mm (\$269). Handcrafted from Premium Grade AISI 304/316 Stainless Steel (Non-Magnetic) blades mounted on a lightweight aluminum guide bar rail, the movable lateral blades measure 17 mm × 22 mm, optimized for precise bone adaptation in standard and delicate thoracic cages. The four movable blades (two per shank) lock at optimum retraction via an integrated turning thread screw lock, while

tapered aluminum handle shanks feature smooth outer edges to protect adjacent soft tissues.

Forged with an anti-glare satin matte finish, Retractor Maker manufactures the Morse-Favaloro Rib Spreader to withstand repetitive hospital autoclave cycles, providing open-heart surgery centers, thoracic clinics, and cardiac surgical theaters worldwide with exceptional mechanical stability, corrosion resistance, and complete autoclave durability.



THANK YOU

We look forward to discussing your surgical retractor requirements.

CONTACT RETRACTOR MAKER

EMAIL

info@retractormaker.com

WHATSAPP

+92 327 8345101

WEBSITES

retractormaker.com
retractormanufacturer.com

ADDRESS

Retractor Maker, Near FFRV+F29, Roras Road, Sialkot-51310, Pakistan

FOLLOW US



[LinkedIn](#)



[Facebook](#)



[Instagram](#)



[X \(Twitter\)](#)



[YouTube](#)



[Pinterest](#)