



RETRACTOR MAKER

PRODUCT BROCHURE OF

Zalkind Malleable Ribbon Retractor

CERTIFIED & APPROVED

✓ EN ISO 13485:2016 + A11:2021



ISO 9001:2015



EN ISO 7153-1:2016



cGMP Approved




SCCI Registered since 2023



SIMAP Approved since 2023

All procedures and processes follow FDA guidelines.



Zalkind Malleable Ribbon Retractor

SKU RM-ZMRR-1

\$44 USD

AVAILABLE SIZES & PRICING

SIZE / OPTION	SKU / DETAILS	PRICE (USD)
10 cm	RM-ZMRR-1-1	\$44
18 cm	RM-ZMRR-1-2	\$49
20 cm	RM-ZMRR-1-3	\$54
21 cm	RM-ZMRR-1-4	\$59
26 cm	RM-ZMRR-1-5	\$64
30 cm	RM-ZMRR-1-6	\$69
33 cm	RM-ZMRR-1-7	\$69

PRODUCT DESCRIPTION

The Zalkind Malleable Ribbon Retractor is an indispensable flexible surgical instrument engineered specifically to provide safe, custom-contoured tissue displacement and visceral shielding across general, plastic, pediatric, and deep orthopedic procedures. Designed with an ultra-ductile flat blade chassis, this spatula allows operating surgeons to shape and mold the instrument intraoperatively to match unique anatomical cavities, avoiding localized pressure trauma or shear injury on delicate visceral organs and neurovascular bundles.

Built to withstand high-frequency autoclave sterilization while offering a lifetime warranty against rust, this Zalkind Malleable Ribbon Retractor series is highly customizable through an extensive matrix of dimensions. Handcrafted from Premium Grade AISI 304/316 Stainless Steel (Non-Magnetic), available overall lengths include 10 cm (\$44), 18 cm (\$49), 20 cm (\$54), 21 cm (\$59), 26 cm (\$64), 30 cm (\$69), and 33 cm (\$69). The blade

widths encompass micro-corridor sizes of 2 mm, 3 mm, 4 mm, 5 mm, 6 mm, and 7 mm, with overview parameters highlighting custom selections stretching up to 2 cm to 7 cm wide for macro tissue displacement. Its smoothly polished edges prevent visceral abrasions while the blade retains its shaped posture under continuous counter-traction.

Forged with an anti-glare satin matte finish, Retractor Maker manufactures the Zalkind Malleable Ribbon Retractor to withstand repetitive manual re-bending and hospital sterilization runs, providing operating rooms, trauma centers, and specialty reconstructive clinics worldwide with exceptional ductility, corrosion resistance, and complete autoclave durability.



THANK YOU

We look forward to discussing your surgical retractor requirements.

CONTACT RETRACTOR MAKER

EMAIL

info@retractormaker.com

WHATSAPP

+92 327 8345101

WEBSITES

retractormaker.com
retractormanufacturer.com

ADDRESS

Retractor Maker, Near FFRV+F29, Roras Road, Sialkot-51310, Pakistan

FOLLOW US



[LinkedIn](#)



[Facebook](#)



[Instagram](#)



[X \(Twitter\)](#)



[YouTube](#)



[Pinterest](#)