



RETRACTOR MAKER

PRODUCT BROCHURE OF

Crawford Aortic Retractor

CERTIFIED & APPROVED

✓ EN ISO 13485:2016 + A11:2021



ISO 9001:2015



EN ISO 7153-1:2016



cGMP Approved




SCCI Registered since 2023



SIMAP Approved since 2023

All procedures and processes follow FDA guidelines.



Crawford Aortic Retractor

SKU RM-CAR-VASC-1

\$44 USD

AVAILABLE SIZES & PRICING

SIZE / OPTION	SKU / DETAILS	PRICE (USD)
Small Model — 19 cm length, 15.0 mm blade	RM-CAR-VASC-1-1	\$44
Large Model — 19 cm length, 15.0 mm blade	RM-CAR-VASC-1-2	\$49

PRODUCT DESCRIPTION

The Crawford Aortic Retractor is an essential cardiovascular surgical retractor engineered specifically to provide safe, non-traumatic displacement of the aorta, vena cava, and major peripheral vessels during vascular reconstruction. Designed with an anatomically curved blade that conforms around large arterial trunks, this instrument allows vascular teams to maintain unhindered operative access while minimizing shear forces on fragile, atherosclerotic vascular walls.

Built for micro-precision in deep thoracic and abdominal fields, this Crawford Aortic Retractor model features a uniform overall length of 19 cm paired with a 15.0 mm blade width. Surgeons can select between two dedicated blade options: a Small Blade (\$44) ideal for procedures requiring high micro-precision within confined surgical fields, and a Large Blade (\$49) configured to elevate broader vascular segments. The balanced ergonomic handle ensures steady, fatigue-free manual counter-traction throughout complex anastomoses.

Manufactured from Premium Grade AISI 304/316 Stainless Steel (Non-Magnetic) with an anti-glare satin matte finish, Retractor Maker produces the Crawford Aortic Retractor to withstand repetitive hospital sterilization cycles, providing cardiovascular, thoracic, and vascular surgery teams with high corrosion resistance, unyielding structural rigidity, and complete clinical durability.



THANK YOU

We look forward to discussing your surgical retractor requirements.

CONTACT RETRACTOR MAKER

EMAIL

info@retractor maker.com

WHATSAPP

+92 327 8345101

WEBSITES

retractor maker.com
retractor manufacturer.com

ADDRESS

Retractor Maker, Near FFRV+F29, Roras Road, Sialkot-51310, Pakistan

FOLLOW US



[LinkedIn](#)



[Facebook](#)



[Instagram](#)



[X \(Twitter\)](#)



[YouTube](#)



[Pinterest](#)