



RETRACTOR MAKER

PRODUCT BROCHURE OF

Heiss Wound Retractor

CERTIFIED & APPROVED

✓ EN ISO 13485:2016 + A11:2021



ISO 9001:2015



EN ISO 7153-1:2016



cGMP Approved



SCCI Registered since 2023



SIMAP Approved since 2023

All procedures and processes follow FDA guidelines.



Heiss Wound Retractor

SKU RM-HWR-1

\$94 USD

AVAILABLE SIZES & PRICING

SIZE / OPTION	SKU / DETAILS	PRICE (USD)
5×5 Prongs Pattern — 10×10 mm or 15×10 mm (Sharp / Blunt)	RM-HWR-1-1	\$94
6×6 Prongs Pattern — 10×10 mm or 15×10 mm (Sharp / Blunt)	RM-HWR-1-2	\$99

PRODUCT DESCRIPTION

The Heiss Wound Retractor is an indispensable minor surgical retractor engineered specifically to provide reliable, hands-free wound edge separation and stable tissue exposure during superficial surgical interventions. By anchoring opposing tissue margins with finely calibrated multi-prong blades, this instrument frees surgical assistants while maintaining a clear, static visual field in delicate operative sites.

Designed for micro-anatomical access, the Heiss Wound Retractor maintains a compact 10.8 cm length framework while offering extensive customizability at checkout. Surgeons can configure the retractor with either a 5×5 prongs pattern (\$94) or a dense 6×6 prongs pattern (\$99), selectable with sharp tips for firm fibrous tissue purchase or blunt tips to prevent localized soft-tissue trauma. Available in 10.0 mm × 10.0 mm and 15.0 mm × 10.0 mm working depths, its precision mechanical thumb-screw allows smooth, micro-incremental arm spreading that stays locked under counter-traction without accidental slipping.

Forged from Premium Grade AISI 304/316 Stainless Steel (Non-Magnetic) with a non-reflective satin matte finish, Retractor Maker manufactures the Heiss Wound Retractor to withstand repeated hospital autoclaving cycles, providing plastic, hand, and general surgery clinics with long-lasting clinical precision.



THANK YOU

We look forward to discussing your surgical retractor requirements.

CONTACT RETRACTOR MAKER

EMAIL

info@retractormaker.com

WHATSAPP

+92 327 8345101

WEBSITES

retractormaker.com
retractormanufacturer.com

ADDRESS

Retractor Maker, Near FFRV+F29, Roras Road, Sialkot-51310, Pakistan

FOLLOW US



[LinkedIn](#)



[Facebook](#)



[Instagram](#)



[X \(Twitter\)](#)



[YouTube](#)



[Pinterest](#)