



RETRACTOR MAKER

PRODUCT BROCHURE OF

Toennis Nerve Root Retractor

CERTIFIED & APPROVED

✓ EN ISO 13485:2016 + A11:2021



ISO 9001:2015



EN ISO 7153-1:2016



cGMP Approved



SCCI Registered since 2023



SIMAP Approved since 2023

All procedures and processes follow FDA guidelines.



Toennis Nerve Root Retractor

SKU RM-TNRR-15-1

\$29 USD

AVAILABLE SIZES & PRICING

SIZE / OPTION	SKU / DETAILS	PRICE (USD)
Toennis Nerve Root Retractor	15 cm overall length, slender shaft, smooth contoured hook tip	\$29
Protective Sleeve / Box	Individual transit packaging	\$29

PRODUCT DESCRIPTION

The Toennis Nerve Root Retractor is a specialized neurosurgical and spinal retractor engineered specifically to provide safe, non-traumatic mobilization and protection of delicate nerve roots and the thecal sac during posterior lumbar discectomy. Designed with an ultra-slender body shaft, this retractor allows spine surgeons to manipulate neural elements gently in narrow interlaminar windows without obstructing microscopic visualization or crowding delicate micro-instruments.

Built to withstand high-frequency autoclave sterilization cycles without degrading, this Toennis Nerve Root Retractor configuration features an overall tool length of 15 cm (\$29) that provides exceptional manual control and refined tactile feedback. Handcrafted from Premium Grade AISI 304/316 Stainless Steel (Non-Magnetic), its working blade profile incorporates a single-ended, slender body layout terminating in a smooth hook contoured specifically for atraumatic neural and delicate tissue manipulation. The ergonomic flat handle ensures steady, fatigue-free manual counter-traction throughout extensive canal decompressions.

Forged with an anti-glare satin matte finish, Retractor Maker manufactures the Toennis Nerve Root Retractor to uphold strict clinical standards, providing neurosurgical suites, spine surgery centers, and orthopedic operating theaters worldwide with exceptional corrosion resistance, structural rigidity, and complete autoclave durability.



THANK YOU

We look forward to discussing your surgical retractor requirements.

CONTACT RETRACTOR MAKER

EMAIL

info@retractormaker.com

WHATSAPP

+92 327 8345101

WEBSITES

retractormaker.com
retractormanufacturer.com

ADDRESS

Retractor Maker, Near FFRV+F29, Roras Road, Sialkot-51310, Pakistan

FOLLOW US



[LinkedIn](#)



[Facebook](#)



[Instagram](#)



[X \(Twitter\)](#)



[YouTube](#)



[Pinterest](#)