



RETRACTOR MAKER

PRODUCT BROCHURE OF

Bookwalter GYN Retractor System

CERTIFIED & APPROVED

✓ EN ISO 13485:2016 + A11:2021



ISO 9001:2015



EN ISO 7153-1:2016



cGMP Approved



SCCI Registered since 2023




SIMAP Approved since 2023

All procedures and processes follow FDA guidelines.

Donald Cameron

Bookwalter GYN Retractor System

SKU RM-BGRS-1

\$1,699 USD


COMPONENT LIST

QUANTITY	COMPONENT	DIMENSIONS / DESCRIPTION
1	Field Post	12" (30 cm)
1	Post Coupler	Standard
1	Horizontal Bar	Standard
4	Tilt Ratchet Mechanisms	Multi-Angle
1	Perineal Ring with Notches	9½" × 14½" (24.1 × 36.8 cm)
2	Vaginal Lateral Blades	1" × 3" (2.5 × 7.6 cm)
1	Vaginal Posterior Blade	1" × 4½" (2.5 × 11.4 cm)
1	Deaver Blade	1" × 5" (2.5 × 12.7 cm)
1	Deaver Blade	1" × 3½" (2.5 × 8.9 cm)
1	Sterilization Box	Protective Case

PRODUCT DESCRIPTION

The Bookwalter GYN Retractor System is a table-mounted self-retaining retractor designed to provide stable, hands-free exposure during gynecological and vaginal procedures. Its modular framework combines a field post, horizontal support bar, notched perineal ring and adjustable tilt ratchet mechanisms to hold selected blades around the operative field.

The GYN configuration includes a 9½" × 14½" notched perineal ring, four multi-angle tilt ratchet mechanisms, vaginal lateral and posterior blades, and two Deaver blades. The complete set is supplied in a dedicated sterilization box for organized storage and processing.

The Bookwalter GYN Retractor System is made from Premium Grade AISI 304/316 Stainless Steel (Non-Magnetic) and is supplied non-sterile. It is designed to be reusable and autoclavable and should be cleaned and sterilized before use according to the healthcare facility's validated procedures.

Retractor Maker supplies the Bookwalter GYN Retractor System directly to hospitals, distributors and other B2B customers.



THANK YOU

We look forward to discussing your surgical retractor requirements.

CONTACT RETRACTOR MAKER

EMAIL

info@retractormaker.com

WHATSAPP

+92 327 8345101

WEBSITES

retractormaker.com
retractormanufacturer.com

ADDRESS

Retractor Maker, Near FFRV+F29, Roras Road, Sialkot-51310, Pakistan

FOLLOW US



[LinkedIn](#)



[Facebook](#)



[Instagram](#)



[X \(Twitter\)](#)



[YouTube](#)



[Pinterest](#)