



RETRACTOR MAKER

PRODUCT BROCHURE OF

Condyle Retractor

CERTIFIED & APPROVED

✓ EN ISO 13485:2016 + A11:2021



ISO 9001:2015



EN ISO 7153-1:2016



cGMP Approved




SCCI Registered since 2023



SIMAP Approved since 2023

All procedures and processes follow FDA guidelines.



Condyle Retractor

SKU RM-CONDR-1

\$44 USD

AVAILABLE SIZES & PRICING

SIZE / OPTION	SKU / DETAILS	PRICE (USD)
8 mm × 27 mm	RM-CONDR-1-1	\$44
10 mm × 44 mm	RM-CONDR-1-2	\$49

PRODUCT DESCRIPTION

The Condyle Retractor is an indispensable maxillofacial surgical retractor engineered specifically to provide safe, slip-resistant exposure of the mandibular condyle and subcondylar neck during orthognathic and trauma procedures. Designed to hook securely behind the posterior border of the ascending ramus, this instrument allows surgeons to displace thick pterygomasseteric tissues while shielding critical neurovascular structures.

Built for deep intraoral access, this Condyle Retractor model features an overall length of 21 cm that provides maximum reach while significantly reducing assistant strain during lengthy surgeries. The instrument is manufactured in two selectable blade dimensions: 8 mm × 27 mm (\$44) and 10 mm × 44 mm (\$49). Its working end incorporates specialized cross-serrations explicitly engineered to provide an enhanced grip on slippery periosteal tissues, ensuring consistent counter-traction without accidental slipping off cortical bone margins.

Machined from Premium Grade AISI 304/316 Stainless Steel (Non-Magnetic) with a reflection-free satin matte finish, Retractor Maker produces the Condyle Retractor to withstand repetitive hospital sterilization cycles, providing oral, maxillofacial, and plastic surgeons worldwide with superior ergonomic balance, high corrosion resistance, and complete autoclave durability.



THANK YOU

We look forward to discussing your surgical retractor requirements.

CONTACT RETRACTOR MAKER

EMAIL

info@retractor maker.com

WHATSAPP

+92 327 8345101

WEBSITES

retractor maker.com
retractor manufacturer.com

ADDRESS

Retractor Maker, Near FFRV+F29, Roras Road, Sialkot-51310, Pakistan

FOLLOW US



[LinkedIn](#)



[Facebook](#)



[Instagram](#)



[X \(Twitter\)](#)



[YouTube](#)



[Pinterest](#)