



RETRACTOR MAKER

PRODUCT BROCHURE OF

Cushing Decompression Retractor

CERTIFIED & APPROVED

✓ EN ISO 13485:2016 + A11:2021



ISO 9001:2015



EN ISO 7153-1:2016



cGMP Approved




SCCI Registered since 2023



SIMAP Approved since 2023

All procedures and processes follow FDA guidelines.



Cushing Decompression Retractor

SKU RM-CDR-215-1

\$59 USD

AVAILABLE SIZES & PRICING

SIZE / OPTION	SKU / DETAILS	PRICE (USD)
Cushing Decompression Retractor	21.5 cm length, 13 × 19 mm fenestrated saddle blade	\$59
Protective Sleeve / Box	Individual transit packaging	\$59

PRODUCT DESCRIPTION

The Cushing Decompression Retractor is an internationally acclaimed neurosurgical and spinal retractor engineered specifically to provide precise, low-pressure tissue retraction and optimal visualization in deep operative corridors. Designed with a distinctive hollow saddle hook geometry, this instrument allows surgeons to mobilize and shield delicate nerve roots, thecal sac margins, and dural tissue while preserving an open visual sightline directly through the blade architecture.

Built to maximize surgical sightlines during intricate decompressions, this Cushing Decompression Retractor model features an overall length of 21.5 cm (\$59) that easily reaches deep spinal and skull base pathology. Its working end incorporates a 13 mm wide by 19 mm long fenestrated blade that secures tissue safely while allowing light and microscopic vision to pass unhindered. The ergonomic flat handle provides surgical assistants with fatigue-free manual leverage and refined tactile feedback.

Forged from Premium Grade AISI 304/316 Stainless Steel (Non-Magnetic) with a reflection-free satin matte finish, Retractor Maker manufactures the Cushing Decompression Retractor to strict clinical standards, delivering high corrosion resistance, unyielding structural rigidity, and complete autoclave compatibility for neurosurgical and spine centers worldwide.



THANK YOU

We look forward to discussing your surgical retractor requirements.

CONTACT RETRACTOR MAKER

EMAIL

info@retractormaker.com

WHATSAPP

+92 327 8345101

WEBSITES

retractormaker.com
retractormanufacturer.com

ADDRESS

Retractor Maker, Near FFRV+F29, Roras Road, Sialkot-51310, Pakistan

FOLLOW US



[LinkedIn](#)



[Facebook](#)



[Instagram](#)



[X \(Twitter\)](#)



[YouTube](#)



[Pinterest](#)