



RETRACTOR MAKER

PRODUCT BROCHURE OF

Wide Self-Retaining Retractor

CERTIFIED & APPROVED

✓ EN ISO 13485:2016 + A11:2021



ISO 9001:2015



EN ISO 7153-1:2016



cGMP Approved



SCCI Registered since 2023



SIMAP Approved since 2023

All procedures and processes follow FDA guidelines.



Wide Self-Retaining Retractor

SKU RM-WSRR-35-1

\$139 USD

AVAILABLE SIZES & PRICING

SIZE / OPTION	SKU / DETAILS	PRICE (USD)
Wide Self-Retaining Retractor	3.5" (8.89 cm) wide smooth blades, 10×12 mm tips, ratchet chassis	\$139
Protective Sleeve / Box	Individual transit packaging	\$139

PRODUCT DESCRIPTION

The Wide Self-Retaining Retractor is a premier self-retaining surgical spreader engineered specifically to provide broad, continuous wound margin displacement and soft-tissue retraction without creating focal pressure points. Designed with expansive solid blades rather than separated prongs, this instrument prevents adipose tissue and muscle bellies from herniating or creeping between teeth, ensuring an unobstructed, static surgical corridor throughout complex procedures.

Built to premium clinical-grade standards (such as New Med Instruments), this Wide Self-Retaining Retractor model features an expansive blade width of 3.5 inches (8.89 cm, \$139) with wide, smooth blades rather than standard separated prongs. Handcrafted from high-quality Premium Grade AISI 304/316 Stainless Steel (Non-Magnetic) to guarantee durability, corrosion resistance, and complete autoclavability, the sleek, low-profile blades measure 10 mm × 12 mm at the working tips to minimize entry trauma while maximizing exposure. The mechanics utilize a positive ratchet locking mechanism combined with ergonomic finger ring handles for smooth, steady, hands-free retraction.

Forged with an anti-glare satin matte finish, Retractor Maker manufactures the Wide Self-Retaining Retractor to withstand repetitive hospital autoclave cycles, delivering unyielding mechanical stability, high corrosion resistance, and complete clinical durability for general surgery operating rooms, orthopedic joint suites, and plastic reconstructive clinics worldwide.



THANK YOU

We look forward to discussing your surgical retractor requirements.

CONTACT RETRACTOR MAKER

EMAIL

info@retractormaker.com

WHATSAPP

+92 327 8345101

WEBSITES

retractormaker.com
retractormanufacturer.com

ADDRESS

Retractor Maker, Near FFRV+F29, Roras Road, Sialkot-51310, Pakistan

FOLLOW US



[LinkedIn](#)



[Facebook](#)



[Instagram](#)



[X \(Twitter\)](#)



[YouTube](#)



[Pinterest](#)