



RETRACTOR MAKER

PRODUCT BROCHURE OF

Cobra Retractor with Hand Rest

CERTIFIED & APPROVED

✓ EN ISO 13485:2016 + A11:2021



ISO 9001:2015



EN ISO 7153-1:2016



cGMP Approved



SCCI Registered since 2023



SIMAP Approved since 2023

All procedures and processes follow FDA guidelines.



Cobra Retractor with Hand Rest

SKU RM-CRHR-1

\$169 USD

AVAILABLE SIZES & PRICING

SIZE / OPTION	SKU / DETAILS	PRICE (USD)
Cobra Retractor with Hand Rest	10.22" length, 32 mm blade, integrated hand rest	\$169
Protective Sleeve / Box	Individual transit packaging	\$169

PRODUCT DESCRIPTION

The Cobra Retractor with Hand Rest is an advanced orthopedic surgical retractor engineered specifically to provide superior mechanical leverage and minimize assistant fatigue during total hip arthroplasty and reconstructive joint surgery. Designed with an integrated ergonomic hand rest platform along its handle, this instrument allows surgical assistants to stabilize their grip and brace their hand, sustaining high counter-traction loads with ease.

Built for smooth anatomical clearance, this Cobra Retractor with Hand Rest model features an overall length of 10.22 inches (\$169) paired with a 32 mm wide blade. Its curved pattern architecture is engineered to facilitate smooth, deep-tissue retraction around femoral and pelvic bone contours, slipping beneath heavy muscular layers without tearing periosteum. The distal tip hooks securely over cortical edges to prevent accidental displacement during vigorous joint reaming.

Forged from premium high-tensile Premium Grade AISI 304/316 Stainless Steel (Non-Magnetic) with a reflection-free satin matte finish, Retractor Maker produces the Cobra Retractor with Hand Rest to withstand intense leverage forces and repeated hospital autoclaving cycles, providing orthopedic surgery teams with unmatched comfort and clinical reliability.



THANK YOU

We look forward to discussing your surgical retractor requirements.

CONTACT RETRACTOR MAKER

EMAIL

info@retractor maker.com

WHATSAPP

+92 327 8345101

WEBSITES

retractor maker.com
retractor manufacturer.com

ADDRESS

Retractor Maker, Near FFRV+F29, Roras Road, Sialkot-51310, Pakistan

FOLLOW US



[LinkedIn](#)



[Facebook](#)



[Instagram](#)



[X \(Twitter\)](#)



[YouTube](#)



[Pinterest](#)