



RETRACTOR MAKER

PRODUCT BROCHURE OF

Willauer-Deaver Retractor

CERTIFIED & APPROVED

✓ EN ISO 13485:2016 + A11:2021



ISO 9001:2015



EN ISO 7153-1:2016



cGMP Approved



SCCI Registered since 2023



SIMAP Approved since 2023

All procedures and processes follow FDA guidelines.



Willauer-Deaver Retractor

SKU RM-WDR-395-1

\$99 USD

AVAILABLE SIZES & PRICING

SIZE / OPTION	SKU / DETAILS	PRICE (USD)
25 mm Wide Blade — 25 mm wide x 108 mm working surface	RM-WDR-395-1-1	\$99
50 mm Wide Blade — 50 mm wide x 108 mm working surface	RM-WDR-395-1-2	\$149

PRODUCT DESCRIPTION

The Willauer-Deaver Retractor is a specialized deep-cavity surgical retractor engineered specifically to provide safe, broad-surface tissue displacement of lung parenchyma, abdominal viscera, and thoracic structures during major open procedures. Designed with an elongated curved hook blade terminating in a rounded non-traumatic lip, this instrument allows surgical assistants to apply sustained counter-traction across heavy visceral planes without causing contusions or serosal tearing.

Built to withstand frequent autoclave sterilization cycles while providing a lifetime warranty against rust, this Willauer-Deaver Retractor configuration features an extraordinary overall tool length of 39.5 cm that provides deep reach into the true pelvis and thoracic cavity. Handcrafted from Premium Grade AISI 304/316 Stainless Steel (Non-Magnetic), the active working blade length measures 108 mm along its curved working surface. The instrument is exceptionally adaptable to varying clinical depths through two distinct choices: a 25 mm wide blade layout (\$99) for focused tissue holding and deep vascular exposure, and an expansive 50 mm wide blade layout (\$149) for broad pulmonary and hepatic displacement. Its ergonomic flat handle provides fatigue-free manual leverage.

Forged with an anti-glare satin matte finish, Retractor Maker manufactures the Willauer-Deaver Retractor to withstand repetitive hospital sterilization cycles, providing cardiothoracic suites, trauma centers, and general surgery operating theaters worldwide with exceptional durability, corrosion resistance, and complete autoclave

durability.



THANK YOU

We look forward to discussing your surgical retractor requirements.

CONTACT RETRACTOR MAKER

EMAIL

info@retractor maker.com

WHATSAPP

+92 327 8345101

WEBSITES

retractor maker.com
retractor manufacturer.com

ADDRESS

Retractor Maker, Near FFRV+F29, Roras Road, Sialkot-51310, Pakistan

FOLLOW US



[LinkedIn](#)



[Facebook](#)



[Instagram](#)



[X \(Twitter\)](#)



[YouTube](#)



[Pinterest](#)