



RETRACTOR MAKER

PRODUCT BROCHURE OF

Biggs Face Lift Retractor, With Fiber optic

CERTIFIED & APPROVED

✓ EN ISO 13485:2016 + A11:2021



ISO 9001:2015



EN ISO 7153-1:2016



cGMP Approved



SCCI Registered since 2023



SIMAP Approved since 2023

All procedures and processes follow FDA guidelines.



Biggs Face Lift Retractor, With Fiber optic

SKU RM-BFLR-FO-1

\$99 USD

COMPONENT LIST

QUANTITY	COMPONENT	DIMENSIONS / DESCRIPTION
1	Biggs Face Lift Retractor, With Fiber optic	9" (229 mm) length, 52 mm blade, fiber-optic bundle
1	Universal Light Cable Adapter	Compatible with Storz, Olympus, Wolf, ACMI fittings
1	Protective Sleeve / Box	Standard transit packaging

PRODUCT DESCRIPTION

The Biggs Face Lift Retractor, With Fiber optic is a premium illuminated plastic surgery instrument engineered to provide exceptional tissue elevation coupled with direct, shadowless lighting during facial reconstructive procedures. By channeling cool illumination directly under undermined skin flaps, this instrument eliminates dark recesses and allows surgeons to clearly identify facial nerve branches and delicate vascular networks.

Designed for precision aesthetic surgery, this Biggs Face Lift Retractor, With Fiber optic model features an overall length of 9 inches (229 mm) and a broad 52 mm wide blade. The blade is contoured with a smooth retaining lip to maintain secure traction on cutaneous tissues without focal pressure necrosis. It integrates an optical fiber bundle featuring a universal connector compatible with standard hospital light cables (Storz, Olympus, Wolf, and ACMI), allowing seamless connection to existing OR illumination sources.

Crafted entirely from Premium Grade AISI 304/316 Stainless Steel (Non-Magnetic) with high-transmission optical fibers and an anti-glare satin finish, Retractor Maker supplies the Biggs Face Lift Retractor, With Fiber optic non-sterile and fully autoclavable, delivering exceptional durability and optical clarity for high-end aesthetic surgery clinics worldwide.



THANK YOU

We look forward to discussing your surgical retractor requirements.

CONTACT RETRACTOR MAKER

EMAIL

info@retractor maker.com

WHATSAPP

+92 327 8345101

WEBSITES

retractor maker.com
retractor manufacturer.com

ADDRESS

Retractor Maker, Near FFRV+F29, Roras Road, Sialkot-51310, Pakistan

FOLLOW US



[LinkedIn](#)



[Facebook](#)



[Instagram](#)



[X \(Twitter\)](#)



[YouTube](#)



[Pinterest](#)