



RETRACTOR MAKER

PRODUCT BROCHURE OF

D'Errico Nerve Root Retractor

CERTIFIED & APPROVED

✓ EN ISO 13485:2016 + A11:2021



ISO 9001:2015



EN ISO 7153-1:2016



cGMP Approved



SCCI Registered since 2023



SIMAP Approved since 2023

All procedures and processes follow FDA guidelines.



D'Errico Nerve Root Retractor

SKU RM-DNRR-1

\$49 USD

AVAILABLE SIZES & PRICING

SIZE / OPTION	SKU / DETAILS	PRICE (USD)
22.8 cm (approx. 9.0") — 7.0 mm wide x 11.0 mm long — Standard lumbar discectomy and nerve root mobilization	RM-DNRR-1-1	\$49
28.0 cm (approx. 11.0") — 7.0 mm wide x 11.0 mm long — Deep operative corridor and bariatric lumbar spine access	RM-DNRR-1-2	\$54

PRODUCT DESCRIPTION

The D'Errico Nerve Root Retractor is a specialized spinal surgical retractor engineered specifically to provide safe, non-traumatic retraction and protection of the thecal sac and traversing nerve roots during posterior lumbar decompressions. Designed with an elongated, slender shaft, this instrument allows spine surgeons to work comfortably in deep interlaminar spaces while keeping manual handles outside the operative line of sight.

Built for delicate neural handling, this D'Errico Nerve Root Retractor model features a fixed 7.0 mm wide by 11.0 mm long angled blade profile. The polished, curved blade tip cradles nerve fibers smoothly, preventing them from slipping back into the operative window or sustaining injury from high-speed burs and rongeurs. Surgeons can select between two reach lengths: 22.8 cm (\$49) for standard spinal corridors or an extended 28 cm (\$54) reach designed for deep or bariatric patient anatomies.

Manufactured from Premium Grade AISI 304/316 Stainless Steel (Non-Magnetic) with a reflection-free satin matte finish, Retractor Maker produces the D'Errico Nerve Root Retractor to strict clinical tolerances, ensuring high corrosion resistance, unyielding structural rigidity, and complete autoclave durability for neurosurgical and spine centers worldwide.



THANK YOU

We look forward to discussing your surgical retractor requirements.

CONTACT RETRACTOR MAKER

EMAIL

info@retractormaker.com

WHATSAPP

+92 327 8345101

WEBSITES

retractormaker.com
retractormanufacturer.com

ADDRESS

Retractor Maker, Near FFRV+F29, Roras Road, Sialkot-51310, Pakistan

FOLLOW US



[LinkedIn](#)



[Facebook](#)



[Instagram](#)



[X \(Twitter\)](#)



[YouTube](#)



[Pinterest](#)