



RETRACTOR MAKER

PRODUCT BROCHURE OF

Cronin Lip & Nostril Retractor

CERTIFIED & APPROVED

✓ EN ISO 13485:2016 + A11:2021



ISO 9001:2015



EN ISO 7153-1:2016



cGMP Approved



SCCI Registered since 2023



SIMAP Approved since 2023

All procedures and processes follow FDA guidelines.

Donald Cameron



Cronin Lip & Nostril Retractor

SKU RM-CLNR-1

\$49 USD

AVAILABLE SIZES & PRICING

SIZE / OPTION	SKU / DETAILS	PRICE (USD)
Cronin Lip & Nostril Retractor	6" (15.2 cm) overall length, 14 mm width, double-ended	\$49
Protective Sleeve / Box	Individual transit packaging	\$49

PRODUCT DESCRIPTION

The Cronin Lip & Nostril Retractor is an indispensable plastic and reconstructive surgery retractor engineered specifically to provide dual-purpose tissue exposure during cleft lip repair and rhinoplasty. Designed with two specialized anatomical working ends on a single chassis, this instrument allows surgeons to transition effortlessly between delicate intranasal expansion and upper labial elevation without switching tools.

Built to maximize intraoperative visualization across two distinct facial structures, this Cronin Lip & Nostril Retractor configuration features an overall length of 6 inches (15.2 cm) and a 14 mm width (\$49). The Nostril End is formatted with two blunt prongs engineered to seat safely within the nasal vestibule, retracting the alar margin to expand the internal visual field. The opposing Lip End features a smooth inward curvature designed to comfortably elevate the upper lip and philtral columns during cleft repair without pinching mucosal tissues.

Forged from Premium Grade AISI 304/316 Stainless Steel (Non-Magnetic) with a reflection-free satin matte finish, Retractor Maker manufactures the Cronin Lip & Nostril Retractor to strict clinical standards, delivering high corrosion resistance, excellent ergonomic balance, and complete autoclave durability for pediatric, plastic, and craniofacial surgery suites worldwide.



THANK YOU

We look forward to discussing your surgical retractor requirements.

CONTACT RETRACTOR MAKER

EMAIL

info@retractormaker.com

WHATSAPP

+92 327 8345101

WEBSITES

retractormaker.com
retractormanufacturer.com

ADDRESS

Retractor Maker, Near FFRV+F29, Roras Road, Sialkot-51310, Pakistan

FOLLOW US



[LinkedIn](#)



[Facebook](#)



[Instagram](#)



[X \(Twitter\)](#)



[YouTube](#)



[Pinterest](#)