



RETRACTOR MAKER

PRODUCT BROCHURE OF

Jackson Vaginal Retractor

CERTIFIED & APPROVED

✓ EN ISO 13485:2016 + A11:2021



ISO 9001:2015



EN ISO 7153-1:2016



cGMP Approved



SCCI Registered since 2023



SIMAP Approved since 2023

All procedures and processes follow FDA guidelines.



Jackson Vaginal Retractor

SKU RM-JVR-1

\$44 USD

AVAILABLE SIZES & PRICING

SIZE / OPTION	SKU / DETAILS	PRICE (USD)
75 mm × 38 mm	RM-JVR-1-1	\$44
90 mm × 38 mm	RM-JVR-1-2	\$49
100 mm × 38 mm	RM-JVR-1-3	\$54

PRODUCT DESCRIPTION

The Jackson Vaginal Retractor is an essential gynecological surgical instrument engineered specifically to provide safe, broad-surface tissue displacement of the vaginal walls and pelvic floor. Designed with an anatomically curved concave blade, this retractor allows surgeons to depress the posterior vaginal wall smoothly, establishing an expansive visual corridor to the cervix and vault without causing focal mucosal tearing.

Built for clinical adaptability, this Jackson Vaginal Retractor model maintains a fixed 18 cm overall length framework while offering multiple selectable model sizes and dynamic blade footprint surface areas from the configuration menu. Surgeons can choose across three standardized blade dimensions: 75 mm × 38 mm (\$44), 90 mm × 38 mm (\$49), and an extended 100 mm × 38 mm (\$54) profile categorized across Size #1 and Size #2 configurations. Its ergonomically angled drop handle provides operating assistants with steady, fatigue-free manual leverage throughout extended pelvic reconstructive procedures.

Forged from Premium Grade AISI 304/316 Stainless Steel (Non-Magnetic) with a reflection-free satin matte finish, Retractor Maker manufactures the Jackson Vaginal Retractor to withstand repetitive hospital sterilization cycles, providing obstetrics, gynecology, and urogynecology teams worldwide with exceptional durability, corrosion resistance, and complete autoclave durability.



THANK YOU

We look forward to discussing your surgical retractor requirements.

CONTACT RETRACTOR MAKER

EMAIL

info@retractormaker.com

WHATSAPP

+92 327 8345101

WEBSITES

retractormaker.com
retractormanufacturer.com

ADDRESS

Retractor Maker, Near FFRV+F29, Roras Road, Sialkot-51310, Pakistan

FOLLOW US



[LinkedIn](#)



[Facebook](#)



[Instagram](#)



[X \(Twitter\)](#)



[YouTube](#)



[Pinterest](#)